1.14402.0001

MQuant®

# Chromate Test

Cr

for the determination of chromium(VI)

## 1. Method

### Determination with color-card comparator

In weakly phosphoric solution chromium(VI) ions react with diphenylcarbazide to form chromium(III) and diphenylcarbazone, which form a red-violet complex. The chromate concentration is measured **semiquantitatively** by visual comparison of the color of the measurement solution with the color fields of a color card.

## 2. Measuring range and number of determinations

| Measuring range / color-scale graduation <sup>1)</sup>                                                               | Number of determinations |
|----------------------------------------------------------------------------------------------------------------------|--------------------------|
| <b>0.005</b> - 0.01 - 0.02 - 0.03 - 0.04 - 0.05 - 0.06 - 0.08 - <b>0.10</b><br>mg/l Cr                               | 150                      |
| <b>0.011</b> - 0.022 - 0.045 - 0.07 - 0.09 - 0.11 - 0.13 - 0.18 - <b>0.22</b><br>mg/l CrO <sub>4</sub> <sup>2-</sup> |                          |

<sup>1)</sup> for conversion factors see section 8

## 3. Applications

This test measures chromium(VI) present in the sample as chromate or dichromate ions.

### Sample material:

Groundwater, surface water, and seawater  
Drinking water  
Industrial water  
Electroplating baths and bath solutions, electroplating wastewater  
Wastewater and percolating water

## 4. Influence of foreign substances

This was checked individually in solutions containing 0.05 and 0 mg/l Cr. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

| Concentrations of foreign substances in mg/l or % |           |                              |           |                                 |             |
|---------------------------------------------------|-----------|------------------------------|-----------|---------------------------------|-------------|
| Al <sup>3+</sup>                                  | 1000      | Fe <sup>3+</sup>             | <b>10</b> | Pb <sup>2+</sup>                | <b>10</b>   |
| Ca <sup>2+</sup>                                  | 50        | Hg <sup>2+</sup>             | <b>10</b> | PO <sub>4</sub> <sup>3-</sup>   | 1000        |
| Cd <sup>2+</sup>                                  | 1000      | Mg <sup>2+</sup>             | 1000      | SiO <sub>3</sub> <sup>2-</sup>  | 1000        |
| CN <sup>-</sup>                                   | 100       | Mn <sup>2+</sup>             | 1000      | Zn <sup>2+</sup>                | 1000        |
| Cr <sup>3+</sup>                                  | 100       | NH <sub>4</sub> <sup>+</sup> | 1000      |                                 |             |
| Cu <sup>2+</sup>                                  | <b>10</b> | Ni <sup>2+</sup>             | 1000      |                                 |             |
| F <sup>-</sup>                                    | 1000      | NO <sub>2</sub> <sup>-</sup> | <b>10</b> |                                 |             |
|                                                   |           |                              |           | EDTA                            | <b>0.1%</b> |
|                                                   |           |                              |           | Surfactants <sup>1)</sup>       | 1%          |
|                                                   |           |                              |           | Na-acetate                      | <b>0.1%</b> |
|                                                   |           |                              |           | NaCl                            | 10%         |
|                                                   |           |                              |           | NaNO <sub>3</sub>               | 10%         |
|                                                   |           |                              |           | Na <sub>2</sub> SO <sub>4</sub> | 10%         |

<sup>1)</sup> tested with nonionic, cationic, and anionic surfactants  
Reducing agents interfere with the determination.

## 5. Reagents and auxiliaries

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

### Package contents:

1 bottle of reagent Cr-1  
3 bottles of reagent Cr-2  
2 test tubes with screw caps (in comparator block)  
1 color card

### Other reagents and accessories:

MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535  
Sodium hydroxide solution 1 mol/l Titripur®, Cat. No. 1.09137  
Sulfuric acid 0.5 mol/l Titripur®, Cat. No. 1.09072  
Chromate standard solution Certipur®, 1000 mg/l CrO<sub>4</sub><sup>2-</sup>, Cat. No. 1.19780  
MQuant® Flat-bottomed long tubes with screw caps for MQuant® tests with color-card comparator (12 pcs), Cat. No. 1.14901

## 6. Preparation

- Analyze immediately after sampling.
- The pH must be within the range 1 - 9.**  
Adjust, if necessary, with sodium hydroxide solution or sulfuric acid.
- Filter turbid samples.

## 7. Procedure

Open the box and set up with both test tubes **on the left**.  
Unfold the color card and insert it, colored end first, into the slit at the lower **right-hand** edge of the box.

|                               | Measurement sample tube <u>nearer to</u> the tester (A)  | Blank tube <u>farther from</u> the tester (B) |                                                                                                   |
|-------------------------------|----------------------------------------------------------|-----------------------------------------------|---------------------------------------------------------------------------------------------------|
| Reagent Cr-1                  | 2 level grey microspoons (in the cap of the Cr-1 bottle) | -                                             | Place into the dry test tube.                                                                     |
| Reagent Cr-2                  | 12 drops <sup>1)</sup>                                   | -                                             | Add, close the tube, and shake <b>vigorously until the reagent Cr-1 is completely dissolved</b> . |
| Pretreated sample (5 - 40 °C) | 20 ml                                                    | 20 ml                                         | Fill the tube to the mark (= 20 ml), close, and mix.                                              |

Leave to stand for 1 min (reaction time).

Slide the color card through to the left until the closest possible color match is achieved between the two open test tubes when viewed from above.  
Read off the result in mg/l Cr or CrO<sub>4</sub><sup>2-</sup> from the color card at the lower right-hand edge of the box.

<sup>1)</sup> Hold the bottle vertically while adding the reagent!

### Notes on the measurement:

- The color of the measurement solution remains stable for at least 60 min after the end of the reaction time stated above.
- Turbidity in the measurement solution makes the color comparison more difficult.
- If the color of the measurement solution is equal to or more intense than the darkest color on the scale, repeat the measurement using **fresh**, diluted samples until a value of less than 0.10 mg/l Cr is obtained.  
Concerning the result of the analysis, the dilution must be taken into account:

$$\text{Result of analysis} = \text{measurement value} \times \text{dilution factor}$$

## 8. Conversions

| Units required                      | = | units given                         | x | conversion factor |
|-------------------------------------|---|-------------------------------------|---|-------------------|
| mg/l CrO <sub>4</sub> <sup>2-</sup> |   | <b>mg/l Cr</b>                      |   | 2.23              |
| <b>mg/l Cr</b>                      |   | mg/l CrO <sub>4</sub> <sup>2-</sup> |   | 0.448             |

## 9. Method control

To check test reagents, measurement device, and handling:  
Dilute the chromate standard solution with distilled water to 0.11 mg/l CrO<sub>4</sub><sup>2-</sup> (0.05 mg/l Cr) and analyze as described in section 7.  
Additional notes see under **www.qa-test-kits.com**.

## 10. Notes

- Reclose the reagent bottles immediately after use.
- Rinse the test tubes **with distilled water only**.
- Information on disposal can be obtained at [www.disposal-test-kits.com](http://www.disposal-test-kits.com).**

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