1.14429.0001

MQuant®

Cyanide Test

CN⁻

for the determination of free cyanide

1. Method

Determination with color-card comparator

Cyanide ions react with a chlorinating agent to form cyanogen chloride, which in turn reacts with 1,3-dimethylbarbituric acid in the presence of pyridine to form a violet dye (König reaction). The cyanide concentration is measured **semiquantitatively** by visual comparison of the color of the measurement solution with the color fields of a color card.

2. Measuring range and number of determinations

Measuring range / color-scale graduation	Number of determinations
0.03 - 0.06 - 0.10 - 0.15 - 0.2 - 0.3 - 0.4 - 0.5 - 0.7 mg/l CN ⁻	200

3. Applications

This test measures only cyanide ions (free cyanide).

Sample material:

Groundwater and surface water
Drinking water and mineral water
Waters from aquaculture
Industrial water
Wastewater and percolating water
Electroplating wastewater
This test is **not suited** for seawater.

4. Influence of foreign substances

This was checked individually in solutions containing 0.3 and 0 mg/l CN⁻. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %					
Ag ⁺	1	F ⁻	1000	Ni ²⁺	0.1
Al ³⁺	100	Fe ³⁺	10	NO ₂ ⁻	10
Ca ²⁺	1000	Hg ²⁺	0.1	Pb ²⁺	100
Cd ²⁺	1000	Mg ²⁺	1000	PO ₄ ³⁻	1000
Cr ³⁺	10	Mn ²⁺	1000	SCN ⁻	0.01
Cu ²⁺	10	NH ₄ ⁺	1000	Zn ²⁺	1000
EDTA	1%	Surfactants ¹⁾	5%	NaCl	1%
		NaNO ₃	10%	Na ₂ SO ₄	10%

All oxidizing and reducing agents interfere with the determination.

¹⁾ tested with nonionic, cationic, and anionic surfactants

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

1 bottle of reagent CN-1
1 bottle of reagent CN-2
1 bottle of reagent CN-3
1 graduated 5-ml plastic syringe
2 test tubes with screw caps (in comparator block)
1 color card

Other reagents and accessories:

MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sodium hydroxide solution 1 mol/l Titripur®, Cat. No. 1.09137
Sulfuric acid 0.5 mol/l Titripur®, Cat. No. 1.09072
MQuant® pH-indicator strips pH 0 - 6.0, Cat. No. 1.09531
Cyanide standard solution Certipur®, 1000 mg/l CN⁻, Cat. No. 1.19533

Refill pack:

Cat. No. 1.18457

Cyanide Test

Refill pack for 1.14798, 1.14429, and 1.14417

(Reagents **without technical accessories** for the number of determinations stated in section 2)

6. Preparation

- Analyze immediately after sampling.
- The pH must be within the range 4.5 - 8.0.**
Adjust, if necessary, with sodium hydroxide solution or sulfuric acid.
- Filter turbid samples.

7. Procedure

Open the box and set up with both test tubes **on the left**.

Slide the comparator block all the way to the left, so that the end holding the test tubes protrudes laterally over the bottom part of the box.

Unfold the color card and insert it, colored end first, into the slit at the lower **right-hand** edge of the box.

	Measurement sample tube <u>nearer</u> to the tester (A)	Blank tube <u>farther</u> from the tester (B)	
Pretreated sample (15 - 40 °C)	5 ml	5 ml	Inject into the test tube with the syringe.
Reagent CN-1	1 level green microspoon (in the cap of the CN-1 bottle)	-	Add, close the tube, and shake vigorously until the reagent is completely dissolved .
Reagent CN-2	1 level green microspoon (in the cap of the CN-2 bottle)	-	Add, close the tube, and shake vigorously until the reagent is completely dissolved . The pH must be within the range 3.0 - 3.5. Check with MQuant® pH-indicator strips.
Reagent CN-3	3 drops ¹⁾	-	Add, close the tube, and mix.

Leave to stand for 5 min (reaction time).

Slide the color card through to the left until the closest possible color match is achieved between the two open test tubes when viewed from above.

Read off the result in mg/l CN⁻ from the color card at the lower right-hand edge of the comparator block within the bottom part of the box.

¹⁾ Hold the bottle vertically while adding the reagent!

Notes on the measurement:

- The color of the measurement solution remains stable for 20 min after the end of the reaction time stated above. (After 30 min the measurement value would have diminished by 3%, after 60 min by 7%.)
- Turbidity in the measurement solution makes the color comparison more difficult.
- If the color of the measurement solution is equal to or more intense than the darkest color on the scale, repeat the measurement using **fresh**, diluted samples until a value of less than 0.7 mg/l CN⁻ is obtained.

Concerning the result of the analysis, the dilution must be taken into account:

$$\text{Result of analysis} = \text{measurement value} \times \text{dilution factor}$$

8. Method control

To check test reagents, measurement device, and handling:

Dilute the cyanide standard solution with distilled water to 0.3 mg/l CN⁻ and analyze as described in section 7.

Additional notes see under **www.qa-test-kits.com**.

9. Notes

- Reclose the reagent bottles immediately after use.
- Rinse the test tubes and the syringe **with distilled water only**.
- Information on disposal can be obtained at www.disposal-test-kits.com.**

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