

MicroCheck™ beverage monitor

Instructions for Use

Introduction

Important

Read these instructions carefully before using the products.

Intended use

The products are intended for research use only, and shall not be used in any clinical or *in vitro* procedures for diagnostic purposes.

Background

MicroCheck™ beverage monitor is a disposable filter unit that is used for recovery and culture of microorganisms from aqueous fluids.

The filter unit has a flexible cylinder that allows for removal of the membrane filter prior to culturing. The product contains a funnel cover that fits onto the funnel base, converting the base into a Petri dish.


The MicroCheck beverage monitor is packaged as two layers of assembled units protected by a box liner bag.

Technical information

Technical data

Support pad:	Cellulose
Cylinder and base:	Polypropylene
Funnel cover:	Polystyrene
Plug:	Polyethylene
Funnel adapter:	Polyethylene
Membrane diameter:	47 mm (1.86 inch)
Filter media:	- ME25 mixed cellulose esters (MCE) membrane, white - Metrical™ Black modified polyether-sulfone (PES) membrane, black
Maximum vacuum:	635 mm Hg (25 in. Hg) (vacuum use only)
Effective filtration area:	13.46 cm ²
Dimensions (H×W):	7.6 × 6.1 cm (3.0 × 2.4 inch)
Sterilization:	Gamma-irradiated

The following illustration shows the parts of the filter unit.




Part	Description
1	Funnel cover
2	Funnel cylinder
3	Membrane
4	Cellulose absorbent pad
5	Funnel base
6	Plug

Instructions for use

Filter the sample

Step	Action
1	Remove the filter unit from the box. Place the filter unit onto the manifold with the applicable adapter assembly as shown below.



Option	Adapter assembly
1	Standard adapter with no. 8 stopper
2	MicroFunnel™ adapter
3	Elongated adapter with no. 8 stopper and auto-clavable adapter

Step	Action
2	Remove the funnel cover without touching the inner surface.
3	(Optional) Rinse the inside walls of the funnel and wet out the membrane surface.
4	Replace the cover and turn on the vacuum pump to draw rinse fluid through. Turn off the vacuum pump.
5	Remove the funnel cover without touching the inner surface and pour the sample into the funnel. Replace the cover.
6	Turn on the vacuum pump for 30 to 40 seconds and draw the sample through the filter.
7	Filter the entire sample first, then rinse the inside walls of the funnel with sterile buffer.
8	Turn off the vacuum pump.

Culture the sample

There are three methods to culture a sample:

- Method A: Remove the membrane and culture the samples on a separate Petri dish.
- Method B: Top load media and culture the samples on the base of the funnel converted into a Petri dish.
- Method C: Bottom load media and culture the samples on the base of the funnel converted into a Petri dish.

Method A

Follow the steps below to remove a membrane and culture the sample on a Petri dish.

Step	Action
1	Hold the filter unit with one hand and gently squeeze the funnel cylinder near the top with the other hand.
2	Remove the funnel cylinder carefully from the base.
3	Dip the tips of a stainless steel forceps into alcohol and flame to sterilize, and allow to cool before handling the membrane filter.
4	Insert the sterilized forceps and remove the membrane from the base leaving behind the support pad.
5	Place the membrane filter on a prepared Petri dish/absorbent pad combination soaked with broth.
6	Transfer the Petri dish to an incubator, and incubate the Petri dish inverted at the specified temperature and time.
7	Count the colonies using your standard procedure.
8	Discard the Petri dish following national and local guidelines for disposal of biohazardous material.

Method B

Follow the steps below to top load media and culture the samples on the base of the funnel converted into a Petri dish.

Step	Action
1	Add media over the surface of the membrane and allow the liquid to be drawn in by residual or momentary vacuum.

Step	Action
2	Remove the filter unit from the manifold, and cap the media port with the supplied plug.
3	Remove the cover and squeeze the funnel gently to separate the cylinder from the base.
4	Disengage the cylinder, and transfer the cover to the base.
5	Discard the funnel cylinder. <i>Result:</i> The base is converted into a Petri dish.
6	Follow steps 6-8 from Method A.

Method C

Follow the steps below to bottom load media and culture the samples on the base of the funnel converted into a Petri dish.

Step	Action
1	Remove the filter unit from the manifold.
2	Slowly add media through the bottom port while rotating the filter unit.
3	Cap the media port with the supplied plug.
4	Remove the cover and squeeze the funnel gently to separate the cylinder from the base.
5	Disengage the cylinder, and transfer the cover to the base.
6	Discard the funnel cylinder. <i>Result:</i> The base is converted into a Petri dish.
7	Follow steps 6-8 from Method A.

Ordering information

For additional ordering information, visit [cytiva.com](https://www.cytiva.com).

Product code	Description	Qty./Pkg.
4761ME	ME25 MCE membrane, gridded, white gridded 0.45 µm 100 mL	50
4763	Metrical Black PES membrane, black gridded, 0.45 µm 100 mL	50
4764	Metrical Black PES membrane, black gridded, 0.8 µm 100 mL	50

Related products

Product code	Description	Qty./Pkg.
4889	Place manifold (quantity 3), including valves (quantity 3), end cap (quantity 1), and hose barb cap (quantity 1)	1 set
4701	Autoclavable funnel adapter	3
82728N	No. 8 rubber stopper	5
4890	MicroFunnel filter funnel adapter	3
4892	Standard adapter	3
4959	Elongated standard adapter	3

Page intentionally left blank



Give feedback on this document

Visit cytiva.com/techdocfeedback or scan the QR code.



cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva.

Metrical, MicroCheck, and MicroFunnel are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Any other third-party trademarks are the property of their respective owners.

© 2024–2025 Cytiva

For local office contact information, visit cytiva.com/contact

29743968 AC V:6 06/2025

