Petrifilm®

6407/6417/6445

Product Instructions

 **(EN)** Yeast and Mold Count Plate

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Product Instructions

Yeast and Mold Count Plate

Product Description and Intended Use

The Neogen® Petrifilm® Yeast and Mold Count (YM) Plate is a sample-ready culture medium system which contains nutrients supplemented with antibiotics, a cold-water-soluble gelling agent, and an indicator that facilitates yeast and mold enumeration. Neogen Petrifilm YM Plates are used for the enumeration of yeast and mold in the food, beverage and bottled water industries. Neogen Petrifilm YM Plate components are decontaminated though not sterilized. Neogen® Petrifilm® Plates are manufactured at an ISO (International Standards Organization) 9001 certified site.

Safety

The user should read, understand, and follow all safety information in the instructions for the Neogen Petrifilm YM Plate. Retain the safety instructions for future reference.

⚠ **WARNING** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage

⚠ WARNING

To reduce the risks associated with exposure to biohazards and environmental contamination:

- Follow current industry standards and local regulations for disposal of biohazardous waste.

To reduce the risks associated with the release of contaminated product:

- Use Neogen Petrifilm YM Plates for food and beverage sample testing that you have validated.
- Follow all product storage instructions contained in the instructions for use.
- Do not use beyond the expiration date.

To reduce risk associated with bacterial infection and workplace contamination:

- Perform Neogen Petrifilm YM Plate testing in a properly equipped laboratory under the control of a skilled microbiologist.
- The user must train its personnel in proper testing techniques. For example, Good Laboratory Practices¹, ISO 7218⁴, or ISO 17025⁵.

To reduce the risks associated with misinterpretation of results:

- Neogen Petrifilm YM Plates do not differentiate any one yeast or mold strain from another.
- Neogen has not documented Neogen Petrifilm YM Plates for use in industries other than food and beverages including bottled water. For example, Neogen has not documented Neogen Petrifilm YM for testing pharmaceuticals, or cosmetics. Neogen has not documented Neogen Petrifilm YM Plates for testing surface and municipal waters, or waters used in the pharmaceutical or cosmetic industries.
- The use of Neogen Petrifilm YM Plates to test water samples in compliance with local water testing regulations is at the sole discretion and responsibility of the end-user. Neogen Petrifilm YM Plates have not been tested with all possible bottled water samples, testing protocols or with all possible strains of microorganisms.
- Do not use Neogen Petrifilm YM Plates in the diagnosis of conditions in humans or animals.

Consult the Safety Data Sheet for additional information.

For information on documentation of product performance, visit our website at www.neogen.com or contact your Neogen representative or authorized distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage

Store unopened Neogen Petrifilm YM Plate pouches refrigerated or frozen at temperatures lower than or equal to 8°C. Just prior to use, allow unopened pouches to come to room temperature before opening (20-25°C / <60% RH). Return unused Neogen Petrifilm YM Plates to pouch. Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool dry place (20-25°C / <60% RH) for no longer than 4 weeks. It is recommended that resealed pouches of Neogen Petrifilm YM Plates be stored in a freezer (see below) if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

To store opened pouches in a freezer, place Neogen Petrifilm YM Plates in a sealable container. To remove frozen Neogen Petrifilm YM Plates for use, open the container, remove the plates that are needed and immediately return remaining plates to the freezer in the sealed container for the remainder of the shelf life. Neogen Petrifilm YM Plates should not be used past their expiration date. The freezer that is used for open pouch storage must not have an automatic defrost cycle as this would repeatedly expose the plates to moisture which can damage the plates.

Do not use plates that show discoloration. Expiration date and lot number are noted on each package of Neogen Petrifilm Plates. The lot number is also noted on individual plates.

⚠ Disposal

After use, Neogen Petrifilm YM Plates may contain microorganisms that may be a potential biohazard. Follow current industry standards for disposal.

For information on potential biohazards, reference Biosafety in Microbiological and Biomedical Laboratories, 6th edition, Section VIII-B: Fungal Agents or equivalent.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Preparation and Incubation of Samples (Bottled Water Excepted)

Sample Preparation

1. Prepare appropriate dilution(s) of the sample as needed.

Use appropriate sterile diluents:

Butterfield's phosphate-buffered dilution water², 0.1% peptone water, peptone salt diluent³, saline solution (0.85-0.90%), Neogen® Wide-Spectrum Neutralizer, bisulphite-free letheen broth or distilled water. **Do not use diluents containing citrate, bisulphite or thiosulfate with Neogen Petrifilm Plates;** they can inhibit growth. If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40-45°C.

2. Blend or homogenize sample.

Plating

1. Place the Neogen Petrifilm YM Plate on a flat, level surface.
2. Lift the top film and dispense 1 mL of sample suspension onto the center of bottom film.
3. Drop the top film down onto the sample.

4. Place the plastic Neogen® Petrifilm® YM Spreader on the center of the plate. Press gently on the center of the spreader to distribute the sample evenly. Spread the inoculum over the entire Neogen Petrifilm YM Plate growth area before the gel is formed. Do not slide the spreader across the film.
5. Remove the spreader and leave the plate undisturbed for at least one minute to permit the gel to form.

Incubation

Incubate Neogen Petrifilm YM Plates for 5 days at 20-25°C in a horizontal position with the clear side up in stacks of no more than 20 plates.

See “Specific Instructions for Validated Methods” for specific requirements.

Preparation and Incubation of Bottled Water Samples

Water filtration and Plate Incubation

1. Following standard procedures for water analysis, membrane filter water sample using a 47 mm, 0.45 micron pore size Mixed Cellulose Ester (MCE) filter.
2. Carefully lift the top film of the Neogen Petrifilm YM Plate. Avoid touching the circular growth area. Place the filter in the center of the plate.
3. Hydrate the Neogen Petrifilm YM Plate by placing 1 mL appropriate sterile, hydration diluent in the center of the filter. Appropriate sterile hydration diluents include distilled water, deionized (DI) water and reverse-osmosis (RO) water.
4. Slowly roll the top film onto the filter. Minimize trapping air bubbles and creating gaps between the filter and the Neogen Petrifilm YM Plate. Lightly apply pressure by using the Neogen Petrifilm Plate YM spreader.
5. Incubate Neogen Petrifilm YM Plates at 20-25°C for 3-5 days⁶ in a horizontal position with the clear side up in stacks of no more than 20.

Interpretation

1. Neogen Petrifilm YM Plates can be counted using a standard colony counter or other illuminated magnifier. Gridlines are visible with the use of a backlight to assist with estimated enumeration.
2. To differentiate yeast and mold colonies on the Neogen Petrifilm YM Plate, look for one or more of the following characteristics:

Yeast	Mold
Small colonies	Large colonies
Colonies have defined edges	Colonies have diffuse edges
Pink-tan to blue-green in color	Variable color
Colonies appear raised (3 dimensional)	Colonies appear flat
Colonies have a uniform color	Colonies have a dark center*

* Mold colonies on the surface of a filter may not exhibit a dark center.

3. Read final yeast and mold results on day 5. Large or fast growing molds may obscure results on Neogen Petrifilm YM Plate by day 5. Check plates on day 3 and record results of plates with high counts (this count can be recorded directly on the plate). If the plate is overgrown by day 5, record the 3-day count as an estimated count.
Note: During the hydration of the Neogen Petrifilm YM Plate with a membrane filter, some colonies may be eluted off of the filter onto the surrounding inoculation area. Count all colonies on both the filter and the surrounding media.
4. Mold colonies may spread and cause the entire growth area to turn blue, black, yellow, etc. Record the three day count as an estimated mold count.
5. High numbers of yeast colonies may cause the entire growth area to turn blue or appear as blue growth around the edge of the inoculated area. If Neogen Petrifilm YM Plates appear to have no growth, lift the top film and examine the gel that adheres to the top film. If numerous yeast are present, you may see white colonies in the gel. This is recorded as a yeast count of too numerous to count (TNTC).
6. The circular growth area is approximately 30 cm². Estimates can be made on plates containing greater than 150 colonies by counting the number of colonies in one or more representative squares and determining the average number per square. Multiply the average number by 30 to determine the estimated count per plate.
7. If a more accurate count is required, re-test the sample plating at higher dilutions.

8. The Neogen Petrifilm YM Plates use a phosphatase enzyme indicator to help detect yeast and mold. All living cells contain phosphatase; therefore natural phosphatase in samples can cause the indicator to react in one of two ways:
 - a) A uniform blue background color (often seen from the organisms used in cultured products).
 - b) Intense, pinpoint blue spots (often seen with spices or granulated products).
9. One or more of the following techniques can help distinguish a color reaction caused by natural phosphatase in a product from yeast and mold colonies:
 - a) Dilute the sample further.
 - b) Allow food particles to settle in the sample, then plate the supernate.
 - c) Check the plate after 24-48 hours of incubation and note any color that is present; if color intensity does not change by day 5 of incubation, the color may be from the phosphatase reaction.
10. Where necessary, colonies may be isolated for further identification. Lift the top film and pick the colony from the gel. If using membrane filtration, the filter may adhere to either the top film or the bottom film when lifting the top film. If the filter adheres to the top film, separate the filter from the top film and pick colonies. Test using standard procedures.
11. If the plates cannot be counted at the end of the 5 day incubation period, store for later enumeration by freezing in a sealable container at temperatures lower than or equal to negative 15°C for no longer than one week.

Note: Delayed counting of Neogen Petrifilm YM Plates with filters is not recommended.

For further information refer to the “Neogen® Petrifilm® Yeast and Mold Count Plate Interpretation Guide.” If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your Neogen Food Safety representative or authorized distributor.

Specific Instructions for Validated Methods

AOAC® Official Methods of AnalysisSM (OMA) #997.02

In an AOAC Official Methods of Analysis study, the Neogen Petrifilm YM Plate method was found to be equivalent to the average log counts of the FDA BAM 7th Edition, Chapter 18 reference method and AOAC SMPR 2021.009.

Scope of Validation:

Enumeration of total yeast and mold in foods and dried cannabis flower.

Incubation:

Incubate Neogen Petrifilm YM Plates for 5 days at 20-25°C.

Interpretation:

Plates containing greater than 150 colonies can either be estimated or recorded as too numerous to count (TNTC). Estimation can be done by counting the number of colonies in one or more representative squares and determining the average number per square. The average number can be multiplied by 30 to determine the estimated count per plate. If a more accurate count is required, the sample can be retested at higher dilutions.

References

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. FDA. Bacteriological Analytical Manual (BAM), Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
3. ISO 6887-1. Microbiology of food and animal feeding stuffs- Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. American Public Health Association. 1998. Standard Methods for the Examination of Water and Wastewater, 20th Ed. Method 9610D.

Refer to the current versions of the standard methods listed above.

Explanation of symbols

info.neogen.com/symbols

AOAC is a registered trademark of AOAC INTERNATIONAL

Official Methods of Analysis is a service mark of AOAC INTERNATIONAL

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
www.neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland

Bray Business Park, Bray, Co. Wicklow,
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2024. All rights reserved.
Neogen and Petrifilm are registered trademarks of Neogen Corporation.
FS00848A

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Références

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. FDA. Bacteriological Analytical Manual (BAM), Reagents Index for BAM consulté sur : <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
3. ISO 6887-1. Microbiology of food and animal feeding stuffs- Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
6. American Public Health Association. 1998. Standard Methods for the Examination of Water and Wastewater, 20th Ed. Method 9610D.

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
Petrifilm[®]

Interpretation Guide

The Neogen[®] Petrifilm[®] Yeast and Mold Count Plate is a sample-ready culture medium system which contains nutrients supplemented with antibiotics, a cold-water-soluble gelling agent, and an indicator that facilitates yeast and mold enumeration. Petrifilm Yeast and Mold Plates are used for the enumeration of yeast and mold in the food and beverage industries.

YM

Yeast and Mold Count Plate

Food and Beverage Applications


Figure 1

Total count = 20
Yeast count = 16
Mold count = 4

Petrifilm Yeast and Mold Count Plate contains both yeast colonies and mold colonies.


Figure 2

Yeast and mold count = 0

Petrifilm Yeast and Mold Count Plate without yeast or molds. Gridlines are visible with the use of a backlight to assist with estimated enumeration.


Figure 3

Estimated total count = 500
Estimated yeast count = 480
Estimated mold count = 21

When colonies number more than 150, estimate the count. Gridlines are visible with the use of a backlight to assist with estimated enumeration. Determine the average number of colonies in one square (1 cm²) and multiply it by 30 to obtain the total count per plate. The inoculated area is approximately 30 cm². Yeast colonies may range in color from tan (as in this example) to pink to blue-green.

For a more accurate count, further dilution of the sample may be necessary.


Figure 4

Estimated yeast count = TNTC

Petrifilm Yeast and Mold Count Plate containing yeast colonies too numerous to count (TNTC). The small, blue colonies at the edge of the plate (highlighted in the box) are present throughout the entire plate although less visible.

For a more accurate count, further dilution of the sample may be necessary.


Figure 5

Estimated mold count = 64

Mold colonies are beginning to crowd and overlap each other on the plate. Count each colony margin or focus. The plate can be divided into sections to assist in counting. In this example, approximately 1/4 of the plate was counted, then the number of colonies counted was multiplied by 4 to get the estimated count on the plate. The section shown has 16 molds.

For a more accurate count, further dilution of the sample may be necessary.


Figure 6a

Mold count = TNTC


Figure 6b

Mold count = 64

Plates in Figures 6a and 6b are the same sample. Figure 6a is a 1:10 dilution and has colonies that are small, faint and numerous, making it difficult to count. Figure 6b is a 1:100 dilution and shows how diluting a sample to obtain a colony count of less than 150 colonies makes counting easier. As with most growth media, in a highly competitive environment (such as Figure 6a), typical colony growth will be inhibited. For heavily contaminated samples such as these, further dilutions are recommended for a more accurate count and more typical colony growth (as in Figure 6b).

Phosphatase Reaction


Figure 7

Yeast and mold count = 0


Figure 8

Yeast and mold count = 0

The Petrifilm Yeast and Mold Count Plates utilize a phosphatase indicator dye. All living cells contain phosphatase; therefore, natural phosphatase in samples can cause the indicator to react. Two types of color reactions are sometimes seen: a uniform blue background color or intense, blue spots. Figure 7 shows uniform blue background color and Figure 8 shows intense blue spots which are often seen with spices or granulated products. Figure 8 also shows food particles that yielded phosphatase.

To reduce a phosphatase reaction, follow one or more of these techniques:

1. Dilute Sample: Further sample dilution will minimize blue background color or reduce the number of intense blue spots.
2. Sample Preparation: Mix sample and let settle before plating. Draw sample from center portion of sample container or use filtered homogenizer bag to avoid plating large particles.
3. Check and Note: Observe plates within 24–48 hours of incubation and make note of any color change to aid in final interpretation.

Bottled Water Applications


Figure 1

Yeast and mold count: 0

Petrifilm Yeast and Mold Count Plate with no colonies.


Figure 2

Yeast and mold count: 10

Mold colonies are large with a dark center and diffuse edge.


Figure 3

Yeast and mold count: 10

Note two small, faint colonies.


Figure 4

Yeast and mold count: 31

Count colonies partially or totally off of the filter.


Figure 5

Yeast and mold count: 51

Estimate colony count when molds merge; darker centers can help enumerate colonies. Count colonies partially or totally off of the filter.

For All Applications

Macroscopic Differentiation


Figure 1

Yeast count = 43

Figure 1 shows typical yeast colonies. Characteristics typical of yeast include:

- Colony is small
- Colony has defined edges
- Colony color can range from pink-tan to blue-green
- Colony may appear raised
- Colony typically is uniform in color, no center focus (dark center)


Figure 2

Mold count = 29

Figure 2 shows typical mold colonies. Characteristics typical of mold include:

- Colony grows large
- Colony has diffuse edges
- Colony color may vary as molds produce a variety of pigments (i.e., brown, beige, orange, blue-green)
- Colony appears flat
- Colony usually has a center focus (i.e., usually darker in color, may also be different color)

Microscopic Differentiation

Yeasts and molds are closely related and cannot always be distinguished from each other without microscopic examination.


Figure 3

To isolate colonies for further identification, lift the top film and pick from the colony within the gel using a loop or similar device.


Figure 4

Transfer the colony to a drop of sterile water on a microscope slide, cover with a coverslip, and view under a microscope.


Figure 5

Yeast typically appear oval and may show budding.


Figure 6

Mold typically appear as branching or thread-like filaments (mycelium).



Figure 7

Molds shown above are in various stages of germination.

Reminders For Use


Food and Beverage Applications

Storage

01


Store the unopened Petrifilm Yeast and Mold Count Plate pouches at refrigerated or frozen at temperatures $\leq 8^{\circ}\text{C}$ ($\leq 46^{\circ}\text{F}$). Use before expiration date on package. Just prior to use, allow unopened pouches to come to room temperature before opening.

02


Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool, dry place for no longer than one month.

Inoculation


03

Place the Petrifilm Yeast and Mold Count Plate on a flat level surface. Lift the top film and with a pipette **perpendicular** to the inoculation area, dispense 1 mL of sample suspension onto the center of the bottom film.


04

Drop the top film down onto the sample.


05

Place the Petrifilm Yeast and Mold Spreader on the center of the plate.

06


Gently apply pressure on the spreader to distribute the inoculum over circular area. **Do not twist or slide the spreader.**

07

Lift the spreader and leave the plate undisturbed for at least one minute to permit the gel to form.


Incubation

08

Incubate plates with clear side up in stacks of up to 20. **Please refer to product instructions for third party validated methods.** Because some molds may grow quickly, it may be useful to read and count plates at 3 days as smaller colonies may be obscured by larger, overgrown molds at 5 days. If this happens, the 3 day count may be used; however, it should be reported as an estimated count.

Interpretation


09

Petrifilm Yeast and Mold Count Plates can be counted using a standard colony counter or other illuminated magnifier.

Reminders For Use


Bottled Water Applications

Storage

01


Store unopened pouches at $\leq 8^{\circ}\text{C}$ ($\leq 46^{\circ}\text{F}$). Use before expiration date on package. Just prior to use, allow unopened pouches to come to room temperature before opening.

02


To seal opened pouches, fold end over and apply adhesive tape. **Do not refrigerate opened pouches.** Use Petrifilm Yeast and Mold Count Plates within one month after opening.

Hydration Procedure


03

Following standard procedures for water analysis, membrane filter water sample using a 47 mm, 0.45 micron pore size Mixed Cellulose Ester (MCE) filter.


04

Place filter in the center of the plate.


05

With the pipette **perpendicular** to the Petrifilm Yeast and Mold Count Plate, place 1 mL of hydration diluent onto the center of the filter. Appropriate sterile diluents include distilled water, deionized (DI) water and reverse osmosis (RO) water.

06


Carefully roll top film down onto the filter.

07

Lightly apply pressure to the Petrifilm Yeast and Mold Spreader to ensure uniform contact of the filter with the gel and to eliminate any air bubbles.


Incubation

08


Incubate Petrifilm Yeast and Mold Count Plates in a horizontal position, clear side up, in stacks on no more than 20 plates at $20\text{--}25^{\circ}\text{C}$ for 3–5 days.

Interpretation

09

Petrifilm Yeast and Mold Count Plates can be counted on a standard colony counter or other illuminated magnifier.

10

Colonies may be isolated for further identification. Lift top film and pick the colony from the gel.

Use Appropriate Sterile Diluents

Butterfield's phosphate buffer, 0.1% peptone water, peptone salt diluent, saline solution (0.85–0.90%), Wide-Spectrum Neutralizer, bisulphite-free letheen broth or distilled water.

Do not use diluents containing citrate, bisulphite or thiosulfate with Petrifilm Yeast and Mold Count Plates; they can inhibit growth.

If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40–45°C.

Neogen offers a full line of products to accomplish a variety of your microbial testing needs.

For more product information, visit info.neogen.com/petrifilm

User's Responsibilities: Neogen Petrifilm Plate performance has not been evaluated with all combinations of microbial flora, incubation conditions and food matrices. It is the user's responsibility to determine that any test methods and results meet the user's requirements. Should re-printing of this Interpretation Guide be necessary, user's print settings may impact picture and color quality.

For detailed CAUTIONS, DISCLAIMER OF WARRANTIES/LIMITED REMEDY and LIMITATION OF NEOGEN LIABILITY, STORAGE AND DISPOSAL information and INSTRUCTIONS FOR USE, see product instructions.

Neogen Corporation, 620 Leshar Place, Lansing, MI 48912 USA.

© Neogen Corporation 2023. All rights reserved. Neogen and Petrifilm are registered trademarks of Neogen Corporation.

FS00582_0823