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Rev.0 / 12.09.2005

Rapid test for the detection of bacterial β -galattosidase

DESCRIPTION

O.N.P.G. TEST is a rapid test for the detection of bacterial β -galattosidase. The tubes contain a desiccated medium suitable for the purpose.

CONTENT OF THE PACKAGES

Each package contains 30 tubes with deiccated medium and one instructions sheet.

ITEMS NECESSARY BUT NOT INCLUDED IN THE PACKAGES

- PHYSIOLOGICAL SOLUTION (code 20095)
- Sundry microbiology laboratory materials

PRINCIPLE OF THE METHOD

The ortho-nitrophenyl-galactopyranoside (O.N.P.G.), contained in the medium, is hydrolyzed by microorganisms able to produce the enzyme β -galactosidase with the formation of a yellow colour orthonitrophenolic compound. Some *Enterobacteriaceae* as for example *E.coli*, *Klebsiella spp.*, *Enterobacter spp.* produce both β -galactosidase and permease and so they are lactose fermentator, others as *Citrobacter spp.* and *Arizona spp.* produce only β -galattosidasi and so they slowly ferment lactose, others as *Salmonella spp.*, *Shigella spp.*, *Proteus spp.*, *Providencia spp.* and *Morganella spp.* do not produce β -galactosidase therefore they are not lactose fermentators. The presence of β -galactosidase enzyme is important for the taxonomy of *Enterobacteriaceae*.

COMPOSITION

Each tube contains a desiccated medium whose formulation includes, among others, the ortho-nitrophenyl-galactopyranoside (O.N.P.G.), which is a specific substratum for the β -galactosidase.

TEST PROCEDURE

- Take the package from the fridge, take one or more tubes and leave it for a few minutes on the bench until it has reached room temperature.
- Add 0.2 mL of PHYSIOLOGICAL SOLUTION (code 20095) to the tubes.
- With a sterile loop suspend in the culture medium of the tube a well isolated bacterial colony of 24 hours incubation chosen from a selective or not selective medium containing lactose.
- Take a second tube and add 0.2 mL of PHYSIOLOGICAL SOLUTION (code 20095). Do not to inoculate this tube and use it as negative control.
- Put the top to the tubes and incubate them at 36 ± 1 °C for a minimum of 4 hours till a maximum of 24 hours.

INTERPRETATION OF THE RESULTS

Interpret the results according to the table 1:

Table 1

Colour of the medium	O.N.P.G TEST
Yellow	Positive
Colourless	Negative

QUALITY CONTROL FOR THE USER

Appearance: amber colour desiccated medium adhering to the tube bottom.

Microbiological control:

Every batch of **O.N.P.G. TEST** is subjected to the quality control, using one bacterial culture of *Citrobacter freundii* ATCC 8090 as positive control and one of *Proteus mirabilis* ATCC 25923 as negative control.

PRECAUTIONS

The product, **O.N.P.G. TEST**, is not classified as hazardous under current legislation, neither does it contain noxious substances in concentrations $\geq 1\%$. **O.N.P.G. TEST** is a disposable device to be used only for diagnostic use *in vitro*. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store **O.N.P.G. TEST** at 2-8 °C in its original package. In such conditions the product is valid until the expiry date shown on the label. Do not use them beyond that date. Dispose of them if they show signs of deterioration.

DISPOSAL OF USED MATERIAL

After the use, **O.N.P.G. TEST** and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.

REFERENCES

1. Edwin H.Lenette: *Manual of Clinical Microbiology* (1995).
2. BLAZEVIC, D.J., and EDERER, G.M.: *Principles of biochemical tests in diagnostic microbiology*. 63-67. New York, John Wiley & Sons, 1975.

PRESENTATION

Product	REF	Σ
O.N.P.G. TEST	88005	30 tubes

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

LIOFILCHEM Bacteriology Products

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745 Fax +390858930330 Website: www.liofilchem.net E-Mail: liofilchem@liofilchem.net

Rev.0 / 12.09.2005