Petrifilm®

6404/6414/6444

## Product Instructions

 (EN) *E. coli* / Coliform Count Plate

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## Product Instructions

### ***E. coli* / Coliform Count Plate**

#### **Product Description and Intended Use**

The Neogen® Petrifilm® *E. coli* / Coliform Count (EC) Plate is a sample-ready-culture medium system which contains modified Violet Red Bile (VRB) nutrients, a cold-water-soluble gelling agent, an indicator of glucuronidase activity, 5-bromo-4-chloro-3-indolyl-D-glucuronide (BCIG), and a tetrazolium indicator that facilitates colony enumeration. Neogen Petrifilm EC Plates are used for the enumeration of *Escherichia coli* (*E. coli*) and coliforms in the food and beverage industries. Neogen Petrifilm EC Plate components are decontaminated though not sterilized. Neogen Food Safety is certified to International Standards Organization (ISO) 9001 for design and manufacturing. Neogen Petrifilm EC Plate has not been evaluated with all possible food products, food processes, testing protocols or with all possible microorganism strains.

#### **Safety**

The user should read, understand, and follow all safety information in the instructions for Neogen Petrifilm EC Plate. Retain the safety instructions for future reference.

**⚠ WARNING:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

#### **⚠ WARNING**

**Do not use this plate for the detection of *E. coli* O157. Because most *E. coli* O157 strains are atypical, for example they are glucuronidase negative, they will not produce a blue color, and will not be detected on Neogen Petrifilm EC Plates.**

#### **To reduce the risks associated with exposure to biohazards and environmental contamination:**

- Follow current industry standards and local regulations for disposal of biohazardous waste.

#### **To reduce the risks associated with release of contaminated product:**

- Follow all product storage instruction contained in the instructions for use.
- Do not use beyond the expiration date.

#### **To reduce the risks associated with bacterial infection and workplace contamination:**

- Perform Neogen Petrifilm EC Plate testing in a properly equipped laboratory under the control of a skilled microbiologist.
- The user must train personnel in current proper testing techniques: for example, Good Laboratory Practices<sup>1</sup>, ISO 17025<sup>2</sup> or ISO 7218<sup>3</sup>.

#### **To reduce the risks associated with misinterpretation of results:**

- Neogen has not documented Neogen Petrifilm EC Plates for use in industries other than food and beverage. For example, Neogen has not documented Neogen Petrifilm EC Plates for testing water, pharmaceuticals, or cosmetics.
- Do not use Neogen Petrifilm EC Plates in the diagnosis of conditions in humans or animals.
- The Neogen Petrifilm EC Plates do not differentiate any one microorganism strain from another.

Consult the Safety Data Sheet for additional information.

If you have questions about specific applications or procedures, please visit our website at [www.neogen.com](http://www.neogen.com) or contact your Neogen representative or authorized distributor.

#### **User Responsibility**

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at [www.neogen.com](http://www.neogen.com), or contact your Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

## Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

## Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

## Storage

Store unopened Neogen Petrifilm EC Plate pouches refrigerated or frozen at temperatures lower than or equal to 8°C (46°F). Just prior to use, allow unopened pouches to come to room temperature before opening. Return unused Neogen Petrifilm EC Plates to pouch. Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool dry place for no longer than four weeks. It is recommended that resealed pouches of Neogen Petrifilm EC Plates be stored in a freezer (see below) if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

To store opened pouches in a freezer, place Neogen Petrifilm EC Plates in a sealable container. To remove frozen Neogen Petrifilm EC Plates for use, open the container, remove the plates that are needed and immediately return remaining plates to the freezer in the sealed container for the remainder of the shelf life. Neogen Petrifilm EC Plates should not be used past their expiration date. The freezer that is used for open pouch storage must not have an automatic defrost cycle as this would repeatedly expose the plates to moisture which can damage the plates.

Do not use Neogen Petrifilm EC Plates that show discoloration. Expiration date and lot number are noted on each package of Neogen Petrifilm EC Plates. The lot number is also noted on individual plates Neogen Petrifilm EC Plates.

## ⚠ Disposal

After use, Neogen Petrifilm EC Plates may contain microorganisms that may be a potential biohazard. Follow current industry standards for disposal.

## Instructions for Use

### Sample Preparation

1. Use appropriate sterile diluents:

Butterfield's phosphate buffered dilution water<sup>4</sup>, 0.1% peptone water, peptone salt diluent, quarter-strength Ringer's solution, saline solution (0.85-0.90%), Neogen® Wide-Spectrum Neutralizer, bisulfite-free letheen broth or distilled water.

**Do not use diluents containing citrate, bisulfite or thiosulfate with Neogen Petrifilm EC Plates; they can inhibit growth.** If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40-45°C.

2. Blend or homogenize sample.
3. For optimal growth and recovery of microorganisms, adjust the pH of the sample suspension to 6.6 - 7.2. For acidic products, adjust the pH with 1N NaOH. For alkaline products, adjust the pH with 1N HCl.

### Plating

1. Place the Neogen Petrifilm EC Plate on a flat, level surface.
2. Lift the top film and with the pipette perpendicular to the inoculation area dispense 1 mL of sample suspension onto the center of bottom film.
3. Roll the top film down onto the sample to prevent trapping air bubbles.
4. Place the Neogen® Petrifilm® Spreader with the flat side down on the center of the plate. Press gently on the center of the Neogen Petrifilm Spreader to distribute the sample evenly. Spread the inoculum over the entire Neogen Petrifilm EC Plate growth area before the gel is formed. Do not slide the Neogen Petrifilm Spreader across the film.
5. Remove the Neogen Petrifilm Spreader and leave the Neogen Petrifilm EC Plate undisturbed for at least one minute to permit the gel to form.

### Incubation

Incubate Neogen Petrifilm EC Plates in a horizontal position with the clear side up in stacks of no more than 20 plates. Several incubation times and temperatures can be used depending on current local reference methods, some of which are listed in the "**Specific Instructions for Validated Methods**" section.

## Interpretation

1. Neogen Petrifilm EC Plates can be counted using a standard colony counter or other illuminated magnifier. Do not count colonies on the foam dam since they are removed from the selective influence of the medium. Do not count artifact bubbles that may be present.

The interpretation of *E. coli* colonies on the Neogen Petrifilm EC Plate is as follows:

AOAC Official Methods (998.08 and 991.14) – enumerate blue to red-blue colonies associated with entrapped gas, regardless of size or intensity of color, as confirmed *E. coli*. Blue colonies without gas are not counted as *E. coli*.

Other coliform colonies are red and closely associated (within one colony diameter) with entrapped gas. Colonies not associated with gas (a distance greater than one colony diameter between colony and gas bubble) are not counted as coliforms. The total coliform count consists of both the red and blue colonies associated with gas at 24 hours. Anytime within the validated method incubation period that a blue colony associated with gas appears, it is a confirmed *E. coli*.

## ▲ WARNING

**Do not use this plate for the detection of *E. coli* O157. Because most *E. coli* O157 strains are atypical, for example they are glucuronidase negative, they will not produce a blue color, and will not be interpreted as *E. coli* on Neogen Petrifilm EC Plates.**

2. The circular growth area is approximately 20 cm<sup>2</sup>. Estimates can be made on plates containing greater than 150 colonies by counting the number of colonies in one or more representative squares and determining the average number per square. Multiply the average number by 20 to determine the estimated count per plate.
3. When present in large numbers, Neogen Petrifilm EC Plates may have one or more of the following characteristics: a deepening of the gel color with many small, indistinct colonies; and many gas bubbles. High concentrations of *E. coli* will cause the growth area to turn blue, while high concentrations of coliforms (non-*E. coli*) will cause the growth area to turn dark red. When this occurs, record results as too numerous to count (TNTC). When an actual count is required, plate at a higher dilution.
4. Where necessary, colonies may be isolated for further identification. Lift the top film using proper testing technique and pick the colony from the gel. Test using standard procedures.
5. If the plates cannot be counted within 1 hour of removal from the incubator, they may be stored for later enumeration by freezing in a sealable container at temperatures lower than or equal to negative 15°C for no longer than one week.

For further information refer to the “Neogen® Petrifilm® *E. coli*/Coliform Count Plate Interpretation Guide”. If you have questions about specific applications or procedures, please visit our website at [www.neogen.com](http://www.neogen.com) or contact your Neogen representative or authorized distributor.

## Specific Instructions for Validated Methods

**AOAC® Official Methods<sup>SM</sup>** (998.08 Confirmed *Escherichia coli* Counts in Poultry, Meats and Seafood, Dry Rehydratable Film Method)

Incubate Neogen Petrifilm EC Plates 24 hours ± 2 hours at 35°C ± 1°C.

**AOAC® Official Methods<sup>SM</sup>** (991.14 Coliform and *Escherichia coli* Counts in Foods, Dry Rehydratable Film Methods)

For coliform results incubate Neogen Petrifilm EC Plates 24 hours ± 2 hours at 35°C ± 1°C.

For *E. coli* results incubate Neogen Petrifilm EC Plates an additional 24 hours ± 2 hours (48 hours ± 4 hours total) at 35°C ± 1°C.

## References

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
4. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
5. ISO 6887. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.

## Explanation of Symbols

[info.neogen.com/symbols](http://info.neogen.com/symbols)

AOAC is a registered trademark of AOAC INTERNATIONAL  
Official Methods is a service mark of AOAC INTERNATIONAL

# Neogen Food Safety

## Neogen Corporation

620 Leshar Place  
Lansing, MI 48912 USA  
Neogen.com

## Neogen Europe Ltd.

The Dairy School  
Auchincruive  
Ayr, KA6 5HU  
Scotland, UK

## Neogen Ireland, Ltd.

Bray Business Park, Bray  
Co. Wicklow  
A98YV29, Ireland

## Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA  
www.neogen.com

© Neogen Corporation 2024. All rights reserved.  
Neogen and Petrifilm are registered trademarks of Neogen Corporation.  
FS00850A

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Petrifilm®

# Interpretation Guide

The Neogen® Petrifilm® *E. coli*/Coliform Count Plate is a sample-ready culture medium system which contains modified Violet Red Bile (VRB) nutrients, a cold-water-soluble gelling agent, an indicator of glucuronidase activity, 5-bromo-4-chloro-3-indolyl-D-glucuronide (BCIG), and a tetrazolium indicator that facilitates colony enumeration. Petrifilm *E. coli*/Coliform Count Plates are used for the enumeration of *Escherichia coli* (*E. coli*) and coliforms in the food and beverage industries.

**EC**

*E. coli*/Coliform Count Plate

The United States Food and Drug Administration (FDA) Bacteriological Analytical Manual (BAM) defines coliforms as Gram-negative rods, which produce acid and gas from lactose fermentation.

Most *E. coli* (about 97%) produce beta-glucuronidase which produces a blue precipitate associated with the colony indicated by the blue to red-blue colonies. The top film traps gas produced by the lactose fermenting coliforms and *E. coli*. About 95% of *E. coli* produce gas, as indicated by colonies associated with entrapped gas (within approximately one colony diameter). Blue colonies without gas are not counted as *E. coli*.<sup>\*</sup> Other coliform colonies are red and closely associated with entrapped gas. The total coliform count consists of both the red and blue colonies associated with gas.<sup>\*</sup>

Most *E. coli* O157 strains are atypical, for example they are glucuronidase negative, they will not produce a blue color, and will not be interpreted as *E. coli* on Petrifilm *E. coli*/Coliform Count Plates.

<sup>\*</sup>Validated through AOAC® Official Methods of Analysis<sup>SM</sup> program. The definition of *E. coli* and coliform varies by regional method. In particular, the confirmation of *E. coli* may vary by country.


Figure 1

***E. coli* count = 49 (blue colonies with gas)**  
**Total coliform count = 87 (red and blue colonies with gas)**

All figures have been counted according to AOAC® Official Methods of Analysis<sup>SM</sup> #998.08 and #991.14. The confirmation of *E. coli* may vary by country.


Figure 2

**Total coliform count = 3 (red colonies with gas)**

Food particles are irregularly shaped and are not associated with gas bubbles.


Figure 3

**Estimated total coliform count = 150**

The recommended counting limit on a Petrifilm *E. coli*/Coliform Count Plate is 150 colonies. Bubble patterns may vary. Gas may disrupt the colony so that the colony “outlines” the bubble. See Circles 1 and 2. Artifact bubbles may result from improper inoculation or from trapped air within the sample. They are irregularly shaped and are not associated with a colony. See Circle 3.

*For a more accurate count, further dilution of the sample may be necessary.*


Figure 4

**No growth = 0**

Notice the changes in gel color in Figures 4 through 10. As the *E. coli* or coliform count increases, the color of the gel turns to dark red or purple-blue. Background bubbles are a characteristic of the gel and are not a result of *E. coli* or coliform growth. See Square 1.


Figure 5

***E. coli* count = 13 (blue colonies with gas)**

**Total coliform count = 28 (red and blue colonies with gas)**

Do not count colonies that appear on the foam barrier because they are removed from the selective influence of the medium. See Circle 1.


Figure 6a


Figure 6b

***E. coli* count = 3 (blue colonies with gas)**

Any blue in a colony (blue to red-blue) with gas indicates the presence of *E. coli*. Front lighting may enhance the detection of blue precipitate formed by a colony. Circle 1 shows a red-blue colony counted using back lighting. Circle 2 shows the same colony with front lighting. The blue precipitate is more evident in Circle 2.


Figure 7

**Estimated *E. coli* = 17 (blue colonies with gas)**

**Estimated total coliform count = 150**

The circular growth area is approximately 20cm<sup>2</sup>. Estimates can be made on plates containing greater than 150 colonies by counting the number of colonies in one or more representative squares and determining the average number per square. Multiply the average number by 20 to determine the estimated count per plate.

*For a more accurate count, further dilution of the sample may be necessary.*

# Too Numerous to Count (TNTC)


Figure 8

## Total coliform count = TNTC

Petrifilm *E. coli*/Coliform Count Plates with colonies that are too numerous to count (TNTC) have one or more of the following characteristics: many small colonies, many gas bubbles and a deepening of the gel color from red to purple-blue.

*For a more accurate count, further dilution of the sample may be necessary.*


Figure 9

## Total coliform count = TNTC

High concentrations of *E. coli* may cause the growth area to turn purple-blue.

*For a more accurate count, further dilution of the sample may be necessary.*


Figure 10

## Total coliform count = TNTC

High concentration of coliforms may cause a deepening of the gel color and many small, indistinct colonies.

*For a more accurate count, further dilution of the sample may be necessary.*


Figure 11


## Total coliform count = TNTC

When high numbers of non-coliform organisms such as *Pseudomonas* are present on Petrifilm *E. coli*/Coliform Count Plates, the gel may turn yellow.

*For a more accurate count, further dilution of the sample may be necessary.*

# Reminders For Use

## Storage

**01**


Store the unopened pouches of plates at frozen or refrigerated temperatures  $\leq 8^{\circ}\text{C}$  ( $\leq 46^{\circ}\text{F}$ ). Use before expiration date on package. It is best to allow pouches to reach room temperature before opening.

**02**


Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store sealed pouches in a cool dry place for no longer than four weeks.

## Inoculation


**03**

Place the Petrifilm *E.coli*/Coliform Count Plate on a level surface. Lift the top film and, with the pipette perpendicular to the inoculation area, dispense 1 mL of sample suspension onto the center of the bottom film.


**04**

Roll the top film down onto sample gently to prevent pushing sample off film and to avoid entrapping air bubbles.


**05**

With flat side down, place Petrifilm Spreader on the center of the Petrifilm *E. coli*/Coliform Count Plate.

**06**


Press firmly on the center of the spreader to distribute the sample evenly. Spread the inoculum over the growth area before the gel is formed. Do not slide the spreader across the top film.

**07**

Remove the spreader and leave the plate undisturbed for one minute to permit the gel to form.


## Incubation

**08**

Incubate Petrifilm *E. coli*/Coliform Count Plates with clear side up in stacks of up to 20. It may be necessary to humidify incubator to minimize moisture loss. **Please refer to the product instructions for third party validated methods.**

## Interpretation

**09**

Petrifilm *E. coli*/Coliform Count Plates can be counted using the Petrifilm Plate Reader Advanced, on a standard colony counter or other illuminated magnifier. Colonies may be isolated for further identification. Lift top film and pick the colony from the gel.

## Use Appropriate Sterile Diluents

Butterfield's phosphate buffered dilution water, 0.1% peptone water, peptone salt diluent, quarter-strength Ringer's solution, saline solution (0.85–0.90%), Wide-Spectrum Neutralizer, bisulfite-free letheen broth or distilled water.


For optimal growth and recovery of the microorganisms, adjust the pH of the sample suspension to 6.6–7.2.

**Do not use diluents containing citrate, bisulfite or thiosulfate with the Petrifilm *E. coli*/Coliform Plates, they can inhibit growth.**

If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40–45°C.

# Bubbles

The illustrations below show examples of various bubble patterns associated with gas producing colonies. All should be enumerated.

Neogen offers a full line of products to accomplish a variety of your microbial testing needs.

For more product information, visit [info.neogen.com/petrifilm](https://info.neogen.com/petrifilm)

**User's Responsibilities:** Neogen Petrifilm Plate performance has not been evaluated with all combinations of microbial flora, incubation conditions and food matrices. It is the user's responsibility to determine that any test methods and results meet the user's requirements. Should re-printing of this Interpretation Guide be necessary, user's print settings may impact picture and color quality.

For detailed CAUTIONS, DISCLAIMER OF WARRANTIES/LIMITED REMEDY and LIMITATION OF NEOGEN LIABILITY, STORAGE AND DISPOSAL information and INSTRUCTIONS FOR USE, see product instructions.

Neogen Corporation, 620 Leshar Place, Lansing, MI 48912 USA.

© Neogen Corporation 2023. All rights reserved. Neogen and Petrifilm are registered trademarks of Neogen Corporation.

FS00580\_0823