Product Instructions

-  (EN) *Bacillus cereus* Count Plate
-  (FR) Plaque de numération de *Bacillus cereus*
-  (DE) *Bacillus cereus* Zählplatte
-  (IT) Piastra per il conteggio di *Bacillus cereus*
-  (ES) Placa de recuento de *Bacillus cereus*
-  (NL) *Bacillus cereus*-tellingplaat
-  (SV) Odlingsplatta för antalsbestämning av *Bacillus cereus*
-  (DA) *Bacillus cereus* tælleplade
-  (NO) *Bacillus cereus* dyrkningsfilm
-  (FI) *Bacillus cereus* kasvatusalusta pesäkkeiden laskentaan
-  (PT) Placa de contagem *Bacillus cereus*
-  (EL) Πλακίδιο Καταμέτρησης *Bacillus cereus*
-  (PL) Płytki do oznaczania liczby *Bacillus cereus*
-  (RU) Тест-пластина для подсчета *Bacillus cereus*
-  (TR) *Bacillus cereus* Sayım Plakası
-  (JA) セレウス菌測定用プレート
-  (ZH) 蜡样芽胞杆菌测试片
-  (TH) แผ่นอาหารเลี้ยงเชื้อสำหรับนับจำนวนบาซิลลัส ซีเรียส
-  (KO) 바실러스 세레우스 카운트 플레이트

Product Instructions

Bacillus cereus Count Plate

Product Description and Intended Use

The Neogen® Petrifilm® *Bacillus cereus* Count (BC) Plate is a selective and differential sample-ready-culture-medium system which contains proprietary nutrients, a cold-water-soluble gelling agent, chromogenic indicators and a lecithinase substrate that facilitates colony enumeration.

This medium is used for the enumeration of the *Bacillus cereus* group, also known as the *B. cereus sensu lato* group, without confirmation in the food and beverage industries. This also includes, but is not limited to, the detection of *B. cytotoxicus*. The Neogen Petrifilm BC Plate components are decontaminated though not sterilized. Neogen Food Safety is certified to International Organization for Standardization (ISO) 9001 for design and manufacturing. The Neogen Petrifilm BC Plate has not been evaluated with all possible food products, food processes, testing protocols or with all possible microorganism strains.

Safety

The user should read, understand, and follow all safety information in the instructions for the Neogen Petrifilm BC Plate. Retain the safety instructions for future reference.

⚠ **WARNING:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ WARNING

To reduce the risks associated with exposure to biohazards and environmental contamination:

- After use, Neogen Petrifilm BC Plates may contain microorganisms that may be a potential biohazard. Follow current local, regional, national and industry standards for disposal.

To reduce the risks associated with release of contaminated product:

- Follow all product storage instructions contained in the instructions for use.
- Do not use beyond the use by date.

To reduce the risks associated with bacterial infection and workplace contamination:

- Perform the Neogen Petrifilm BC Plate testing in a properly equipped laboratory under the control of a skilled microbiologist.
- The user must train its personnel in current proper testing techniques: for example, Good Laboratory Practices¹, ISO 7218², or ISO 17025³.

To reduce the risks associated with misinterpretation of results:

- Neogen has not documented the Neogen Petrifilm BC Plates for use in industries other than food and beverage. For example, Neogen has not documented the Neogen Petrifilm BC Plates for testing water, pharmaceuticals, or cosmetics.
- Acceptance of the Neogen Petrifilm BC Plate method for the testing of water per an accepted local government regulation is at the sole discretion and responsibility of the end user.
- Do not use the Neogen Petrifilm BC Plates in the diagnosis of conditions in humans or animals.
- The Neogen Petrifilm BC Plates do not differentiate any one *Bacillus cereus sensu lato* strain from another.

For information on documentation of product performance, visit our website at www.neogen.com or contact your Neogen representative or authorized Neogen distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your Neogen representative or authorized Neogen distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please call Customer Service or your official Neogen Food Safety representative for a Returned Goods Authorization.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. IN NO EVENT SHALL NEOGEN'S LIABILITY UNDER ANY LEGAL THEORY EXCEED THE PURCHASE PRICE OF THE PRODUCT ALLEGED TO BE DEFECTIVE.

Storage

Store the unopened Neogen Petrifilm BC Plate pouches refrigerated or frozen at temperatures (-20 to 8°C/-4 to 46°F). Just prior to use, allow the unopened Neogen Petrifilm BC Plate pouches to come to room temperature before opening (20-25°C (68-77°F)/<60% RH). Return the unused Neogen Petrifilm BC Plates to pouch. Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed Neogen Petrifilm BC Plate pouches in a cool dry place for no longer than four weeks. It is recommended that resealed pouches of Neogen Petrifilm BC Plates be stored in a freezer for no longer than four weeks if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

To store opened pouches in a freezer, place the Neogen Petrifilm BC Plates in a sealable container. To remove the frozen Neogen Petrifilm BC Plates for use, open the container, remove the plates that are needed, and immediately return the remaining plates to the freezer in the sealed container for no longer than four weeks. The freezer that is used for open pouch storage must not have an automatic defrost cycle as this would repeatedly expose the Neogen Petrifilm BC Plates to moisture, which can damage the plates.

Do not use Neogen Petrifilm BC Plates that show discoloration. Use by date and batch code are noted on each package of the Neogen Petrifilm BC Plates. The lot number is also noted on individual Neogen Petrifilm BC Plates. The Neogen Petrifilm BC Plates should not be used past their use by date.

△ Disposal

After use, Neogen Petrifilm BC Plates may contain microorganisms that may be a potential biohazard. Follow current local, regional, national and industry standards for disposal.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Sample Preparation

1. Use appropriate sterile diluents:

Butterfield's phosphate buffer⁴, 0.1% peptone water, peptone salt diluent⁵, phosphate buffered saline, buffered peptone water⁵, saline solution (0.85-0.90%), bisulfite-free letheen broth, or distilled water.

Do not use diluents containing citrate, bisulfite, or thiosulfate with the Neogen Petrifilm BC Plates; they can inhibit growth. See section "Specific Instructions for Validated Methods" for specific requirements.

2. Blend or homogenize the sample.
3. For optimal growth and recovery of microorganisms in acidic products (<pH 5), adjust the pH of the sample suspension to a pH greater than pH 5 using 1N NaOH.

Plating

1. Place the Neogen Petrifilm BC Plate on a flat, level surface.
2. Lift the top film, and with the pipette perpendicular to the inoculation area, dispense 1 mL of sample suspension onto the center of bottom film.
3. Roll the top film down onto the sample to prevent trapping air bubbles.
4. Place the Neogen® Petrifilm® Flat Spreader (6425) with the flat side down on the center of the plate. Press gently on the center of the spreader to distribute the sample evenly. Spread the inoculum over the entire Neogen Petrifilm BC Plate growth area before the gel is formed. Do not slide the spreader across the film.
5. Remove the Neogen Petrifilm Flat Spreader and leave the plate undisturbed for at least one minute to permit the gel to form.

Incubation

Incubate the Neogen Petrifilm BC Plates in a horizontal position with the clear side up in stacks of no more than 10 plates.

Incubate the Neogen Petrifilm BC Plates for 20-24 hours at 30°C±1°C or 35°C±1°C.

Interpretation

1. The Neogen Petrifilm BC Plates can be counted on a benchtop or using a standard colony counter or other illuminated magnifier. Do not count colonies on the foam dam since they are removed from the selective influence of the medium.
2. Count all red-violet colonies with a cream/white precipitate around the colony as *Bacillus cereus*. These typical colonies may or may not be associated with a yellow acid zone.
3. Non-*Bacillus cereus* organisms will be inhibited or may appear as blue colonies or pinpoint red colonies with no cream/white precipitate around the colony.
4. Certain strains of *Bacillus cytotoxicus* may occasionally exhibit a darker red-purple color with a distinct white precipitate at 30°C±1°C with typical morphology and appearance at 35°C±1°C.
5. The circular growth area is approximately 30 cm². Estimates can be made on the Neogen Petrifilm BC Plates containing greater than 100 colonies. Count the number of colonies in one or more representative squares and determine the average number per square. Multiply the average number by 30 to determine the estimated count per Neogen Petrifilm BC Plate.
6. High concentrations of *Bacillus cereus* on the Neogen Petrifilm BC Plates may appear as many small, indistinct colonies and/or darkening of the outer edge. When this occurs, record results as TNTC. For a more accurate count, further dilution of the sample may be necessary.
7. Where necessary, colonies may be isolated for further identification. Lift the top film and pick the colony from the gel. Test using standard procedures.
8. If the Neogen Petrifilm BC Plates cannot be counted within the incubation period, they may be stored for later enumeration by freezing in a sealable container at temperatures lower than or equal to -15°C (5°F) for no longer than one week.

For further information refer to the “Neogen® Petrifilm® *Bacillus cereus* Count Plate Interpretation Guide”. If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your representative or authorized Neogen distributor.

References

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO 7218, Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025, General requirements for the competence of testing and calibration laboratories.
4. FDA Bacteriological Analytical Manual (BAM), Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
5. ISO 6887, Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension, and decimal dilutions for microbiological examination.

6. ISO 7932, Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of presumptive *Bacillus cereus* - Colony-count technique at 30 degrees C
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.

Explanation of Symbols

info.neogen.com/symbols

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
Neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2024. All rights reserved.
Neogen and Petrifilm are registered trademarks of Neogen Corporation.
FS01537A

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Referenser

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO 7218, Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025, General requirements for the competence of testing and calibration laboratories.
4. FDA Bacteriological Analytical Manual (BAM), Reagents Index för BAM finns på:
<http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
5. ISO 6887, Mikrobiologi i livsmedel och djurfoder – Beredning av testprover, initial suspension och decimalspädningar för mikrobiologisk undersökning.
6. ISO 7932, Mikrobiologi i livsmedel och djurfoder - Horisontell metod för att räkna upp presumtiva *Bacillus cereus* - Koloniräkningsteknik vid 30 °C
7. ISO 16140-2. Mikrobiologi i livsmedelskedjan – Metodvalidering – Protokoll för validering av alternativ (proprietära) metoder mot en referensmetod.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO 7218, Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025, General requirements for the competence of testing and calibration laboratories.
4. FDA Bacteriological Analytical Manual (BAM), Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
5. ISO 6887, Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension, and decimal dilutions for microbiological examination.
6. ISO 7932, Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of presumptive *Bacillus cereus* - Colony-count technique at 30 degrees C
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO 7218, Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025, General requirements for the competence of testing and calibration laboratories.
4. FDA Bacteriological Analytical Manual (BAM), BAM용 시약 색인: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
5. ISO 6887, Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension, and decimal dilutions for microbiological examination.
6. ISO 7932, Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of presumptive *Bacillus cereus* - Colony-count technique at 30 degrees C
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.

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
Figure 13

Bacillus cereus count = 75

Food particles are irregularly shaped or filamentous (Circle 1 and 2). Do not enumerate food particles.

Artifact bubbles may result from improper inoculation or from trapped air within the sample. They are irregularly shaped and are not associated with a colony (Circle 3). Do not enumerate artifact bubbles.


Figure 14

Bacillus cereus count = 65

Do not count colonies on the foam barrier (Circle 1) because they are removed from the selective influence of the medium.

Colonies in close proximity of one another may share a cream/white zone (Circle 2) and should be counted as two colonies.


Figure 15

Bacillus cereus count = 49

Bacillus cytotoxicus may occasionally exhibit a darker red-purple color with a distinct cream/white zone at 30°C±1°C (Figure 15) with typical morphology and appearance at 35°C±1°C (Figure 16) .



Figure 16

Bacillus cereus count = 31

Reminders For Use

Storage


01
Store the **unopened** Petrifilm *Bacillus cereus* Count Plate pouches at frozen or refrigerated temperature equal to -20 to 8°C (-4 to 46°F). Use before expiration date on package. It is best to allow pouches to reach room temperature before opening.


02
Seal by folding the end of the pouch over and applying adhesive tape. To prevent exposure to moisture, do not refrigerate opened pouches. Store resealed pouches in a cool dry place (20–25°C/<60% RH) for no longer than four weeks.


Inoculation


03
Place the Petrifilm *Bacillus cereus* Plate on **level** surface. Lift the top film and with the pipette **perpendicular** dispense 1 mL of sample suspension onto the center of bottom film.


04
Roll top film down onto sample **gently** to prevent pushing sample off film and to avoid entrapping air bubbles.

05
Place the Petrifilm Flat Spreader on the center of the Petrifilm *Bacillus cereus* Count Plate.


06
Gently apply pressure on spreader to distribute inoculum over circular area. **Do not** twist or slide the spreader.


07
Lift spreader. Wait a minimum of one minute for gel to solidify.

Incubation

08
Incubate plates with clear sides up in stacks up to 10. **Please refer to the product instructions for third party validated methods.**

Interpretation

09
Petrifilm *Bacillus cereus* Count Plates can be counted on a benchtop, or using a standard colony counter or other illuminated magnifier.

Use Appropriate Sterile Diluents

Butterfield's phosphate buffer, buffered peptone water, 0.1% peptone water, peptone salt diluent, saline solution (0.85–0.90%), bisulfite-free letheen broth or distilled water.

Do not use diluents containing citrate, bisulfite or thiosulfate with the Petrifilm *Bacillus cereus* Count Plates; they can inhibit growth.

If needed, adjust the pH of the sample suspension to a pH greater than pH 5.

User's Responsibilities: Neogen Petrifilm Plate performance has not been evaluated with all combinations of microbial flora, incubation conditions and food matrices. It is the user's responsibility to determine that any test methods and results meet the user's requirements. Should re-printing of this Interpretation Guide be necessary, user's print settings may impact picture and color quality.

Neogen offers a full line of products to accomplish a variety of your microbial testing needs.

For more product information, visit info.neogen.com/petrifilm

Neogen Corporation, 620 Leshar Place, Lansing, MI 48912 USA.

© Neogen Corporation 2024. All rights reserved. Neogen and Petrifilm are registered trademarks of Neogen Corporation.

For detailed CAUTIONS, DISCLAIMER OF WARRANTIES/LIMITED REMEDY and LIMITATION OF NEOGEN LIABILITY, STORAGE AND DISPOSAL information and INSTRUCTIONS FOR USE, see product instructions.


FS01457_0724

Petrifilm®

Interpretation Guide

The Neogen® Petrifilm® *Bacillus cereus* Count (BC) Plate is a selective and differential sample-ready-culture-medium system which contains nutrients, a cold-water-soluble gelling agent, chromogenic indicators and a lecithinase substrate that facilitates colony enumeration of *Bacillus cereus* group, also known as the *Bacillus cereus sensu lato* group, in as soon as 20 hours.


Figure 1

***Bacillus cereus* Count = 0**

Figure 1 shows a *Bacillus cereus* Count Plate without colonies.


Figure 2

***Bacillus cereus* Count = 54**

Count all red-violet colonies with a cream/white zone as *Bacillus cereus*.


Figure 3

***Bacillus cereus* count = 69**

Bacillus cereus colonies may occasionally be associated with a yellow acid zone in addition to the cream/white zone. Count all red-violet colonies with a cream/white zone as *Bacillus cereus*.


Figure 4

***Bacillus cereus* count = 10**

Figure 4 shows a *Bacillus cereus* Count Plate with a few bacterial colonies.

Too Numerous Too Count (TNTC)

Petrifilm *Bacillus cereus* Count Plates that are TNTC may have one or more of the following characteristics: many small indistinct colonies, darkening of the outer edge and deepening of the gel color from blue-purple to dark purple.


Figure 5a

***Bacillus cereus* count = 0**

High concentrations of colonies on the plates will cause the entire growth area to turn blue-purple to deep purple.

For a more accurate count, further dilution of sample may be necessary.


Figure 5b

***Bacillus cereus* count = TNTC**


Figure 7

***Bacillus cereus* count = TNTC**

The counting range on a Petrifilm *Bacillus cereus* Count Plate is less than or equal to 100 colonies.

High concentrations of colonies on the plate may appear as many small indistinct colonies on the plate.

For a more accurate count, further dilution of sample may be necessary.


Figure 6

***Bacillus cereus* count = TNTC**

The counting range on a Petrifilm *Bacillus cereus* Count Plate is less than or equal to 100 colonies.

High concentrations of colonies on the plate may appear as many small indistinct colonies on the plate.

For a more accurate count, further dilution of sample may be necessary.


Figure 8

***Bacillus cereus* count = TNTC
Estimated count = ~ 360**

The circular growth area is approximately 30cm². Gridlines are visible with the use of a backlight to assist with estimated enumeration. Estimates can be made on Petrifilm *Bacillus cereus* Count Plates by counting the number of colonies in two or more representative squares and determining the average number per square. Multiply the average number by 30 to determine the estimated count per plate.

For a more accurate count, further dilution of sample may be necessary.


Figure 9

***Bacillus cereus* count = 8**

Bacillus cereus may appear as large colonies. Count all red-violet colonies with a cream/white zone regardless of size and color intensity.


Figure 10

***Bacillus cereus* count = 29**

Bacillus cereus colonies may be small (Circle 1) or large (Circle 2) in size. Count all red-violet colonies with a cream/white zone regardless of size and color intensity.


Figure 11

***Bacillus cereus* count = 0**

Non-*Bacillus cereus* colonies may be inhibited or may appear as red colonies without a cream/white zone (Figure 11) or as blue colonies (Figure 12). Do not interpret red colonies without a cream/white zone or blue colonies as *Bacillus cereus*.


Figure 12

***Bacillus cereus* count = 0**