





Description	5ml Tips A Type	5ml Tips B Type
Maximum Volume	5000µl	5000µl
Graduation	1000µl, 2000µl, 3000µl, 4000µl, 5000µl	1000µl, 2000µl, 3000µl, 4000µl, 5000µl
Material	Virgin Polypropylene	Virgin Polypropylene
Autoclavable	Yes	Yes
Certified RNase, DNase, DNA and Pyrogen free	Yes	Yes
Hazards Data	Health 0 Non-hazardous Flammability 1 Flammable (Class A fire) Reactivity 0 Non-reactive *Consult the polymer MSDS for further information on potential hazards associated with this material.	
Production Standards & Procedures	Product is Produced under controlled environment. All personnel are capped, masked, gloved and gowned in the assembly area. Lot size is variable due to continuous molding process. This Product is manufactured to the following specifications: Magnified visual examination for flash (+0.02in).	

Packaging	AT-5000A  Bulk Bulk-Sterile	AT-5000B  Bulk Bulk-Sterile
-----------	--	--

Technical Drawing	Tips 5000 A 	Tips 5000 B 
-------------------	--	--

Order in:- AT-5000A, AT-5000B

Volume	Bulk	Sterile	Individually Wrapper
5000µl	B	S	W

B=Bulk, S=Sterile, W=Wrapped