

PEROXYDE 0-100PPM BTE 100 BANDELETTES

Référence

0299003

Low Level Peroxide Test Strip, 100ppm

1. Remove one test strip from the vial. Do not touch the test pad with your fingers.
2. Dip the test strip into the solution being tested for 1 second.
3. Remove and shake off any excess liquid.
4. Compare the test pad to the color chart after 10 seconds.

0 ppm



1 ppm



3 ppm



10 ppm



50 ppm



100 ppm



Ces spécifications ne peuvent être considérées que comme des informations générales et ne dispenseront pas l'utilisateur de son obligation de s'assurer que le produit est adapté pour l'utilisation prévue. En cas de doute, veuillez nous contacter.

FT 0299003

avril 24