PDS No. 6289x0	PRODUCT DATA SHEET			Page 1 of 1	
Dovision 02	Cell Culture Dish, CELLCOAT®, 60 x 15 mm				
Revision 03	Greiner Item-No. 6289x0				greiner bio-one
Valid for Item-No.:	628910	628920	628930	628940	628950

1.	Description / Specification		
1.1	Description	Cell Culture Dish, CELLCOAT®, 60 x 15 mm, with vents	
	·	628910: Laminin protein coating	
		628920: Fibronectin protein coating	
		628930: Poly-L-Lysine protein coating	
		628940: Poly-D-Lysine protein coating	
		628950: Collagen Type I protein coating	
1.2	Dimensions	See Customer Drawing	
		Total weight: 7,4 – 7,8 g	
1.3 Volume Max. volume: 28 ml		Max. volume: 28 ml	
		Working volume: 6-7 ml	
		Growth area: 21 cm <sup>2</sup>	
1.4 Material / Resin <u>Dish</u> : 628910: PS (Polystyrene), coa		Dish: 628910: PS (Polystyrene), coated with Laminin	
		628920: PS (Polystyrene), coated with Fibronectin	
		628930: PS (Polystyrene), coated with Poly-L-Lysine	
		628940: PS (Polystyrene) coated with Poly-D-Lysine	
		628950: PS (Polystyrene), coated with Collagen Type I	
		Dish / Lid: the materials for manufacturing are free of heavy metals	
1.5	Colour	<u>Dish</u> : clear	
		<u>Lid</u> : clear	
1.6	Sterilisation	No	
1.7	7 Quality Control - <u>Raw Material-Control</u> : physical testing		
		- Product-Control: testing of attributive and variable characteristics in	
		accordance with the valid specification	
1.8	Other Information	- For single use only	
		- Elevations for ventilation of culture	

2.	Features	
2.1	Basic features	-
2.2	Temperature range	628910, -920: +2℃ to +8℃
		628930, -940, -950: room temperature
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage:
		https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	628910, -920: 6 month after month of production
		628930, -940, -950: 2 years after month of production
2.7	Other Information	-

3.	Packaging	628910, -920	628930, -940, -950
3.1	Pieces / Bag	5	20
3.2	Pieces / Box	20	100
3.3	Lot-No.	AABB0D89 (year, week, country code, index, unique code)	AABB0D89 (year, week, country code, index, unique code)
3.4	Other Information	-	-

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Bata Chock dubject to charige Without Hotico.					
Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this	
Revision	Date	Date	Date	document or drawing is confidential and proprietory to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.	
02	2 March 2015	3 March 2015	3 March 2015		
Date	Name	Name	Name		
18.02.2013	S. Kaelberer	Dr. T. Schreiber	A. Schulz		