PDS No. 6579x0	PRODUCT DATA SHEET			Page 1 of 1	
Boyisian 06	6 Well Cell Culture Multiwell Plate, CELLCOAT®				9
Revision 06	Greiner Item-No. 6579x0				greiner bio-one
Valid for Item-No.:	657910	657920	657930	657940	657950

1.	Description / Specification	
1.1	Description	PS Plate, 6 well, with vents, solid F-bottom (flat), alphanumeric well coding,
		single position lid_with condensation rings
		657910: Laminin protein coating
		657920: Fibronectin protein coating
		657930: Poly-L-Lysine protein coating
		657940: Poly- D -Lysine protein coating
		657950: Collagen Type I protein coating
1.2	Dimensions	See Customer Drawing
1.3	Total volume: 16,1 ml (mathematically calculated)	
		Working volume: 2 - 5 ml
		Growth area: 9,6 cm ²
1.4	Material / Resin	Plate: 657910: PS (Polystyrene) coated with Laminin
		657920: PS (Polystyrene) coated with Fibronectin
		657930: PS (Polystyrene) coated with Poly-L-Lysine
		657940: PS (Polystyrene) coated with Poly- D -Lysine
		657950: PS (Polystyrene) coated with Collagen Type I
		Lid: PS (Polystyrene)
		The materials for manufacturing are free of heavy metals
1.5	Colour	Plate: clear
		<u>Lid</u> : clear
1.6	Sterilization	No
1.7	.7 Quality Control - Raw Material-Control: physical testing	
		- Product-Control: testing of attributive and variable characteristics in
		accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	657910, -920: +2°C to +8°C 657930, -940, -950: room temperature
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	4800 x g: swinging-bucket rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html (Only concerning the standard plate without coating)
2.6	Shelf life	657910, -920: 6 month after month of production 657930, -940, -950: 2 years after month of production
2.7	Other Information	-

3.	Packaging	657910, -920	657930, -940, -950
3.1	Pieces / Bag	1	5
3.2	Pieces / Box	5	50
3.3	Lot-No.	AABB0D89 (year, week, country code, index, unique code)	AABB0D89 (year, week, country code, index, unique code)
3.4	Other Information	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this	
Revision	Date	Date	Date	document or drawing is confidential and proprietory to Greiner Bio-One GmbH. This	
05	17 March 2017	23 March 2017	24 March 2017	document may not be reproduced for any	
Date	Name	Name	Name	reason without written permission from Greiner Bio-One GmbH. All rights of design, invention,	
02.12.2014	S. Kaelberer	Dr. T. Schreiber	A. Schulz	and copyright are reserved.	