PDS No. 6649x0	PRODUCT DATA SHEET				Page 1 of 1
Dovision 04	Cell Culture Dish, CELLCOAT®, 100 x 20 mm				6
Revision 04	Greiner Item-No. 6649x0				greiner bio-one
Valid for Item-No.:	664910	664920	664940	664950	

1.	Description / Specification			
1.1 Description		Cell Culture Dish, CELLCOAT®, 100 x 20 mm with vents, heavy design 664910: Laminin protein coating		
		664920: Fibronectin protein coating		
		664940: Poly-D-Lysine protein coating		
		664950: Collagen Type I protein coating		
1.2	Dimensions	See Customer Drawing		
		Weight: lid: 8,1 – 8,6 g		
		dish: 10,3 – 10,8 g		
1.3 Volume Max. volume: 100 ml		Max. volume: 100 ml		
		Working volume: 16 - 17 ml		
		Growth area: 58 cm <sup>2</sup>		
1.4	Material / Resin	Dish: 664910: PS (Polystyrene), coated with Laminin		
		664920: PS (Polystyrene), coated with Fibronectin		
		664940: PS (Polystyrene) coated with Poly-D-Lysine		
		664950: PS (Polystyrene), coated with Collagen Type I		
		Dish / Lid: the materials for manufacturing are free of heavy metals		
1.5	Colour	Dish: clear		
		Lid: clear		
1.6	Sterilization	No		
1.7	Quality Control Raw Material-Control: physical testing			
		Product-Control: testing of attributive and variable characteristics in		
		accordance with the valid specification		
1.8	Other Information	- For single use only		
		- Elevations for ventilation of culture		

2.	Features	
2.1	Basic features	-
2.2 Temperature range 664910, -920: +2℃ to +8℃		664910, -920: +2℃ to +8℃
		664940, -950: room temperature
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage:
		https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	664910, -920: 6 month after month of production
		664940, -950: 2 years after month of production
2.7	Other Information	-

3.	Packaging	664910	664920	664940, -950
3.1	Pieces / Bag	10	5	10
3.2	Pieces / Box	10	10	40
3.3	Lot-No.	AABB0D89 (year, week, country code, index, unique code)	AABB0D89 (year, week, country code, index, unique code)	AABB0D89 (year, week, country code, index, unique code)
3.4	Other Information	-		-

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this
Revision	Date	Date	Date	document or drawing is confidential and proprietory to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
03	8 April 2015	8 April 2015	8 April 2015	
Date	Name	Name	Name	
04.03.2015	S. Kaelberer	Dr. T. Schreiber	A. Schulz	