

501002 Deep Well Plates

Features

1. USP VI Polypropylene material
2. High chemical stability;
3. Autoclavable
4. Centrifugal strength up to 3000g
5. Compatible with multi-channel pipette and automatic equipment
6. Optimized bottom design to minimize residual volume
7. Clear alphanumeric codes for easy identification.
8. Conform to SBS/ ANSI
9. Specific 96 well deep plate used in nucleic acid extraction
10. Stackable, reliable tight sealing
11. Non-pyrogenic, DNase/RNase free
12. Both E-beam sterilized and non-sterilized deep well plates are available

Description

Temperature Range : Approx.-80°C to+120°C
 Shelf Life : 5 Years after date of production (ensure package is in good)
 Sterile : NO

Product Range

Cat.No.	Volume (mL)	Well Number	Well Type	Shape Type	/Pack	/Case
501002	1.0	96	Square Well	U	5	50

Technical Drawing

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© Wuxi NEST Biotechnology Co., Ltd. Version 1, Jan. 10th, 2018