

PCR Sealing Films

Nest PCR sealing film can be used with both Elisa plates and PCR plates. The film can minimize the solution evaporation even under high temperature. Both common and advanced films are available.

Features

Common Sealing Film Product Features

1. Non-permeable film, medical grade glue
2. Compatible with most PCR plates
3. Good sealing performance

Advanced PCR Sealing Film Product Features

1. Non-tacky to skin and gloves
2. Minimal to no Auto-fluorescence

Description

Temperature Range : Common Sealing Film : suitable for use between -20°C—110°C
 Advanced PCR Sealing Film : suitable for use between -70°C—100°C

Shelf Life :
 Sterile : No

Product Range

Cat.No.	Name	Spec (mm)	/Pack	/Case
410001	PCR Common Sealing Film	146 × 81	100	500
410011	Advanced sealing film for qPCR	141 × 78	100	500

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA