

Microcentrifuge Tubes

Introduction

Microcentrifuge tubes are made of transparent high polymer materials – polypropylene (PP), used with microcentrifuge and widely used in molecular biology, clinical chemistry and biochemical research.

Features

1. Using imported PP materials, resistant to chemical corrosion and low temperature, compliant to USP Class-6 standard
2. Ultra-high centrifugal stability, centrifugal tolerance up to 30000xg
3. Endotoxin unit less than 0.1 EU
4. Withstand sterilization under high temperature and pressure 121°C/15psi
5. Tube cap can be opened or closed with one hand, easy to operate
6. Available in sterile or unsterile options
7. None of DNA\RNA enzyme or endotoxin

Materials

Tube: PP

Product Range

Microcentrifuge tube, unsterilized

Cat No.	Volume (ml)	Type	Sterile	Piece/Box	Box/Carton
605001	0.6	Transparent, cone bottom	No	1000	10
615001	1.5	Transparent, cone bottom	No	500	10
620011	2	Transparent, round bottom	No	500	10

**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
Email : info@nest-wuxi.com
Website : www.cell-nest.com
Add : No.530, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**


Tel : +1-732-381-0268
Email : info@nestscientificusa.com
Website : www.nestscientificusa.com
Add : 1592 Hart St., Rahway, NJ07065, USA

For Research Use Only. © Wuxi NEST Biotechnology Co., Ltd.

Microcentrifuge tube, sterilized

Cat No.	Volume (ml)	Type	Sterile	Piece/Box	Box/Carton	Cat No.
605601	0.6	Transparent, cone bottom	Yes	50	750	10
615601	1.5	Transparent, cone bottom	Yes	50	400	10
620611	2	Transparent, round bottom	Yes	50	400	10
603011	5	Transparent, round bottom	No	200	1	10
603111	5	Transparent, round bottom	Yes	200	1	10
605401	0.6	Brown, cone bottom	Yes	50	15	10
615401	1.5	Brown, cone bottom	Yes	50	8	10

Technical Drawing



**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
Email : info@nest-wuxi.com
Website : www.cell-nest.com
Add : No.530, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**

Tel : +1-732-381-0268
Email : info@nestscientificusa.com
Website : www.nestscientificusa.com
Add : 1592 Hart St., Rahway, NJ07065, USA