



Specification

1.04410.0050 Mercury(II) acetate for analysis EMSURE® ACS,Reag. Ph Eur

	Spec. Values	
Assay (complexometric)	≥ 99.0	%
Substances insoluble in dilute acetic acid	≤ 0.01	%
Chloride (Cl)	≤ 0.005	%
Nitrate (NO ₃)	≤ 0.005	%
Sulfate (SO ₄)	≤ 0.005	%
other heavy metals (as Pb)	≤ 0.002	%
Fe (Iron)	≤ 0.001	%
Mercury (I) (as Hg)	≤ 0.3	%
Residue after reduction	≤ 0.02	%

Corresponds to ACS,Reag. Ph Eur

Claudia Wiegand

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.