



Specification

1.04874.0250 Potassium hydrogen phthalate for analysis EMSURE® Reag. Ph Eur

	Spec. Values	
Assay (Perchloric acid titration from dried substance)	≥ 99.5	%
Identity (IR-spectrum)	passes test	
In water insoluble matter	≤ 0.005	%
Chloride (Cl)	≤ 0.002	%
Sulfur compounds (as SO ₄)	≤ 0.005	%
Heavy metals (as Pb)	≤ 0.0005	%
Cu (Copper)	≤ 0.0002	%
Fe (Iron)	≤ 0.0005	%
Na (Sodium)	≤ 0.01	%
Pb (Lead)	≤ 0.0005	%
Loss on Drying (105 °C)	≤ 0.2	%

Corresponds to Reag. Ph Eur

Claudia Wiegand

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.