



# Specification

---

1.04984.0500 Potassium hexacyanoferrate(II) trihydrate for analysis EMSURE® ACS,ISO,  
Reag. Ph Eur

---

	Spec. Values	
Assay (cerimetric)	99.0 - 102.0	%
Insoluble matter	≤ 0.005	%
Chloride (Cl)	≤ 0.01	%
Sulfate (SO <sub>4</sub> )	≤ 0.005	%
Na (Sodium)	≤ 0.01	%
Pb (Lead)	≤ 0.002	%

Corresponds to ACS,ISO,Reag. Ph Eur

Claudia Wiegand

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.