PDS No. 112101	PRODUCT DATA SHEET	Page 1 of 1
Davisian 05	Tube, 4 ml, PS, Round Bottom	6
Revision 05	Greiner Item-No. 112101	greiner

1.	Description / Specification			
1.1	Description	PS tube with round bottom, 4 ml		
		Appropriate grip stopper, CatNo. 303321		
1.2	Dimensions	See Customer Drawing		
1.3	Volume	Working volume: 3 ml		
1.4	Material / Resin	PS (Polystyrene), free of heavy metal		
1.5	Colour	Clear		
1.6	Sterilization	No		
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification		
1.8	Intended Use	General laboratory product for the processing, centrifugation and storage of samples to be used by qualified personnel in a laboratory environment.		
1.9	Other Information	For single use only		

2.	Features	
2.1	Basic features	-
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation	5.800 x g: swinging-bucket rotor 6.200 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	240
3.2	Pieces / Box	3.600
3.3	Lot-No.	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this
Revision 04	Date 18 August 2021	Date 23.09.2021	Date 25.10.2021	document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any
Date 17.04.2015	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.