



A CALIBRE SCIENTIFIC COMPANY



Reagecon pH 7.00 Colour Coded Buffer Solution at 20°C

Item No. 10705C



Use Colour Coded pH buffer for ease of use. Reagecon's pH Colour Coded Buffer Solution tested at 20°C is tested by an ISO 17025 Accredited Test Method (A2LA Ref: 6739.03) and NIST traceable; has guaranteed stability throughout its entire shelf life; even after opening the bottle.

Additional Details

Pack Size	500ml
Weight	0.62
Weight Unit	Kg
Colour	Yellow
Tariff Code	3822 190090
Country of Origin	Ireland

DESCRIPTION

NIST traceable

PRODUCT SPECIFICATIONS

pH Buffer Standard, pH Buffer Material

Value	pH 7.00 ± 0.01
Additional Information	Accuracy ± 0.01 pH units
Tested	at 20°C

Use	Laboratory Chemical
Packaging	HDPE Bottle GL45
Physical Form	Liquid