


PDS No. 38108x	PRODUCT DATA SHEET			Page 1 of 1
Revision 04	EVA CapMat (96 square naps)			 greiner bio-one
	Greiner Item-No. 38108x			
Valid for Item-No.:	381080	381081 (sterile)		

1.	Description / Specification	
1.1	Description	EVA CapMat for sealing 96 well plates, square naps, alphanumeric nap coding 381080: for 2 ml 96 well MASTERBLOCK®, not pierceable 381081: for 2 ml 96 well MASTERBLOCK®, not pierceable, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	N/A
1.4	Material / Resin	Ethylene Vinyl Acetate (EVA)
1.5	Colour	Translucent
1.6	Sterilization	381080: no 381081: SAL 10 ⁻³
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	381081: free of detectable DNase/RNase, human DNA and pyrogens 381080, -081: resistant to DMSO
2.2	Temperature range	-60°C to +65°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	Weak acids: none Strong acids: oxidizing acids attack Weak alkalis: none Strong alkalis: none Organic solvents: resistant below 80°C
2.6	Shelf life	381080: n/a 381081: 5 years after month of production
2.7	Other Information	-

3.	Packaging	381080	381081
3.1	Pieces / Bag	10	1
3.2	Pieces / Box	50	50
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-	Certificate of Quality

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 17 April 2015	Date 20 April 2015	Date 20 April 2015	
Date 11.12.2009	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.