

PDS No. 5420x0	<b>PRODUCT DATA SHEET</b>			Page 1 of 1
Revision 04	<b>EASYstrainer for 50 ml Tubes</b>			
	Greiner Item-No. 5420x0			
Valid for Item-No.:	<b>542000</b> (sterile)	<b>542040</b> (sterile)	<b>542070</b> (sterile)	

1.	Description / Specification	
1.1	Description	EASYstrainers are Cellstrainers consist of a slightly conical polypropylene housing, with a filtration tissue in three color-codes mesh sizes - fits all standard 50 ml tubes - handle and mantle surface for improved aseptic handling - vent opening for fast filtration - stackable - sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	N/A
1.4	Material / Resin	Housing: PP (Polypropylene), free of heavy metal Filter: PET, free of heavy metal
1.5	Colour	Housing: 542000: yellow 542040: green 542070: blue Filter: natural
1.6	Sterilization	SAL 10 <sup>-3</sup>
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	Warehouse conditions: +3°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Blister bag	1
3.2	Pieces / Box	50
3.3	Lot-No.	M YY ZZZZ (M=Batch-No., YY=year, ZZZZ= consecutive counter) or E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 22 July 2020	Date 23 July 2020	Date 23 July 2020	
Date 09.07.2019	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.