



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	165KS01
GENERAL DESCRIPTION	WOODEN APP.B/01PCS 1000PKG
DETAILED DESCRIPTION	1 box containing 1000 paper-film pouches each containing 1 wooden applicator with cotton tip.

PRODUCT COMPONENTS

APPLICATOR	WOODEN APPLICATOR	DIMENSIONAL SPECIFICATIONS
		TOTAL LENGTH: 150.50 mm \pm 2.00 mm STICK DIAMETER (*purchased product, no tolerances): 2.30 mm TIP DIAMETER: 4,8 mm \pm 1,0 mm TIP LENGTH: 15.0 mm \pm 2.0 mm STICH LENGTH (*purchased product, no tolerances): 148.0 mm

PACKAGING SPECIFICATION

OUTER BOX QTY	1000
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.01-21/07/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.