



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	167KS01
GENERAL DESCRIPTION	PL.APPL.RAYON IND.WRAP.1000PCS
DETAILED DESCRIPTION	1 box containing 1000 paper-film pouches each containing 1 plastic applicator with rayon tip

PRODUCT COMPONENTS

APPLICATOR	PLASTIC APPLICATOR	DIMENSIONAL SPECIFICATIONS
		STICK LENGTH: 144.0 mm \pm 1.0 mm PLASTIC TIP DIAMETER: 2.50 mm \pm 0.20 mm FIBRE TIP LENGTH: 15.0 mm \pm 2.0 mm FIBRE TIP DIAMETER: 4.50 mm \pm 0.50 mm TOTAL LENGTH: 147.0 mm \pm 1.0 mm

PACKAGING SPECIFICATION

PACKAGING	PAPER-FILM POUCH
BOX QTY	1000

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.01-13/06/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.