



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

| | |
|-----------------------------|---|
| PRODUCT CODE | 947C |
| GENERAL DESCRIPTION | SRK 4ML UNS BULK 100MM 300PKG |
| DETAILED DESCRIPTION | 6 boxes containing 50 tubes with green screw cap containing a plastic applicator with rayon tip and 4ml of medium UNS |

PRODUCT COMPONENTS

| MEDIUM | COMPOSITION |
|---------------|--|
| UNS MLS | Peptone Sodium Chloride Disodium phosphate Monopotassium phosphate Tween 80 Lecithin Sodium thiosulfate Sodium Thioglicollate Sodium Bicarbonate Potassium chloride |

| SAMPLING TUBE | TUBE WITH ROUNDED BOTTOM | DIMENSIONAL SPECIFICATIONS |
|----------------------|---------------------------------|---|
| | | TUBE LENGTH: 101.80 mm \pm 0.20 mm TUBE DIAMETER: 16.10 mm \pm 0.20 mm INTERNAL DIAMETER TUBE: 14.10 mm \pm 0.20 mm |
| CAP | GREEN CAP | DIMENSIONAL SPECIFICATIONS |
| | | CAP DIAMETER: 21.50 mm \pm 0.20 mm CAP LENGTH: 16.30 mm \pm 0.20 mm |

This bulletin (Rev.04-20/05/2015) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

| APPLICATOR | PLASTIC APPLICATOR | DIMENSIONAL SPECIFICATIONS |
|------------|--------------------|--|
| | | STICK LENGTH: 83.30 mm \pm 0.50 mm DIAMETER HANDLE PART: 2.50 mm \pm 0.20 mm TIP LENGTH: 16.50 mm \pm 2.00 mm FIBRE TIP DIAMETER: 4.50 mm \pm 0.50 mm TOTAL LENGTH: 85.50 mm \pm 2.00 mm |

PACKAGING SPECIFICATION

| | |
|---------------|------|
| INNER BOX QTY | 50 |
| OUTER BOX QTY | 300 |
| PACKAGING | BULK |

IMPORTANT INFO

| | |
|----------------------|----------------------|
| STERILIZATION METHOD | GAMMA IRRADIATION |
| STORAGE TEMPERATURE | 5/25°C-DO NOT FREEZE |
| SHELF LIFE | 12 MONTHS |
| MEDIUM pH | 7.3 \pm 0.5 |