



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	948C
GENERAL DESCRIPTION	SRK 4ML UNS BULK 160MM TUBE
DETAILED DESCRIPTION	1 outer box contain 6 inner boxes each composed of 50 tubes with green cap filled with 4mL of UNS solution and 1 rayon tip swab each ones.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UNS MLS	Peptone Sodium Chloride Disodium phosphate Monopotassium phosphate Tween 80 Lecithin Sodium thiosulfate Sodium Thioglicollate Sodium Bicarbonate Potassium chloride Osmotic water

APPLICATOR	PLASTIC APPLICATOR	SPECIFICATIONS
		TOTAL LENGTH: 147.00 mm $\pm$ 2.00 mm STICK DIAMETER: 2.50 mm $\pm$ 0.10 mm FIBER TIP LENGTH: 16.00 mm $\pm$ 3.00 mm FIBER TIP DIAMETER: 4.50 mm $\pm$ 0.50 mm STICK MATERIAL: POLYSTYRENE TIP MATERIAL: RAYON
SAMPLING TUBE	TUBE	SPECIFICATIONS
		TOTAL LENGTH: 158.50 mm $\pm$ 0.20 mm TUBE DIAMETER: 16.00 mm $\pm$ 0.20 mm INTERNAL DIAMETER: 14.50 mm $\pm$ 0.20 mm

This bulletin (Rev.02-17/01/2022) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.

**TECHNICAL DATA SHEET**

CAP	GREEN CAP	SPECIFICATIONS
		LOWER CAP DIAMETER: 21.50 mm $\pm$ 0.20 mm CAP LENGTH: 16.30 mm $\pm$ 0.20 mm HIGHER CAP DIAMETER: 20.50 mm $\pm$ 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

 PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	BULK

 IMPORTANT INFO

STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	12 MONTHS
MEDIUM pH	7.3 $\pm$ 0.5

This bulletin (Rev.02-17/01/2022) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.