

TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	155C
GENERAL DESCRIPTION	TAMPONE SECCO A/PL 1000 PKG
DETAILED DESCRIPTION	10 boxes each containing 100 tubes with white cap with inside a plastic applicator with rayon tip in bulk

PRODUCT COMPONENTS

APPLICATOR	PLASTIC APPLICATOR	DIMENSIONAL SPECIFICATIONS
		STICK LENGTH: 144.0 mm \pm 1.0 mm PLASTIC TIP DIAMETER: 2.50 mm \pm 0.20 mm FIBRE TIP LENGTH: 15.0 mm \pm 2.0 mm FIBRE TIP DIAMETER: 4.50 mm \pm 0.50 mm TOTAL LENGTH: 147.0 mm \pm 1.0 mm
CAP	WHITE CAP	DIMENSIONAL SPECIFICATIONS
		EXTERNAL CAP DIAMETER: 12.00 mm \pm 0.10 mm CAP LENGTH: 32.80 mm \pm 0.20 mm INSERTION DIAMETER: 9.30 mm \pm 0.05 mm
TRANSPORT TUBE	TUBE	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 154.50 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.2 mm EXTERNAL DIAMETER: 12.0 mm \pm 0.5 mm

This bulletin (Rev.02-27/06/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.

**TECHNICAL DATA SHEET****PACKAGING SPECIFICATION**

INNER BOX QTY	100
OUTER BOX QTY	1000
PACKAGING	BULK

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.02-27/06/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.