



TECHNICAL DATA SHEET

1) PRODUCT DESCRIPTION

PRODUCT CODE:	150C
GENERAL DESCRIPTION:	DRY SWAB WOODEN APPLICATOR 1000PKG
DETAILED DESCRIPTION:	10 Boxes each containing 100 tubes with 1 wooden applicator

2) PRODUCT COMPONENTS

APPLICATOR	DIMENSIONAL SPECIFICATIONS	COMPONENTS	RAW MATERIALS
WOODEN APPLICATOR	LENGTH OF THE WOODEN SHAFT: 148 mm \pm 1 mm STICK DIAMETER HANDLE PART: 2.25 mm \pm 0.5 mm LENGTH OF THE FIBRE TIP: 15 mm \pm 3 mm FIBER TIP DIAMETER: 4.80 mm \pm 1.0 mm TOTAL LENGTH: 150.5 mm \pm 2mm	WOODEN STICK	WOOD
		FIBER	COTTON
		GLUE	VINILIC ADHESIVE

TUBE

TUBE COMPONENTS	DIMENSIONAL SPECIFICATIONS	RAW MATERIALS
CAP	EXTERNAL DIAMETER: 12 mm \pm 1 mm CAP LENGTH: 32.5 mm \pm 1 mm	POLYETHILEN HIGH DENSITY
TUBE	EXTERNAL DIAMETER: 12 mm \pm 0.5 mm INTERNAL DIAMETER: 9.5 mm \pm 0.5 mm TUBE LENGTH: 155 mm \pm 1 mm TOTAL LENGTH (CAP+TUBE): 175 mm \pm 2 mm	POLYPROPYLEN

This bulletin is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



3) PACKAGING SPECIFICATION

INNER BOX QTY:	100
OUTER BOX QTY:	1000
PACKAGING:	BULK

4) OTHER SPECIFICATION

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin is for your guidance and is based on information believed to be reliable and should be periodically reviewed.