



## TECHNICAL DATA SHEET

### 1) PRODUCT DESCRIPTION

<b>PRODUCT CODE:</b>	170KS01
<b>GENERAL DESCRIPTION:</b>	ALLUMINIUM APP.B/01PCS 1000PKG
<b>DETAILED DESCRIPTION:</b>	1 Box containing 1000 Individually Paper-Film pouched alluminium applicators

### 2) PRODUCT COMPONENTS

APPLICATOR	DIMENSIONAL SPECIFICATIONS	COMPONENTS	RAW MATERIALS
ALLUMINIUM APPLICATOR	LENGTH OF THE ALLUMINIUM SHAFT: 148 mm $\pm$ 1 mm ALLUMINIUM TIP DIAMETER: 1.0 mm $\pm$ 0.2 STICK DIAMETER HANDLE PART: 1.0 mm $\pm$ 0.2 LENGTH OF THE FIBRE TIP: 15 mm $\pm$ 3 mm FIBER TIP DIAMETER: 2.50 mm $\pm$ 1.0 TOTAL LENGTH: 150 mm $\pm$ 2mm	ALLUMINIUM SHAFT	ALLUMINIUM
		FIBER	RAYON
		GLUE	VINILIC ADHESIVE

### 3) PACKAGING SPECIFICATION

<b>BOX QTY:</b>	1000
<b>PACKAGING:</b>	PAPER-FILM POUCH
<b>TOP SHEET MATERIAL:</b>	MEDICAL GRADE PAPER PLUS 60GR
<b>BOTTOM SHEET MATERIAL:</b>	PA/PE FILM - 100 $\mu$

### 4) OTHER SPECIFICATION

<b>STERILIZATION METHOD</b>	GAMMA IRRADIATION
<b>STORAGE TEMPERATURE</b>	2/30°C – DO NOT FREEZE
<b>SHELF LIFE</b>	36 MONTHS

This bulletin is for your guidance and is based on information believed to be reliable and should be periodically reviewed.