

TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	160C
GENERAL DESCRIPTION	TAMPONE SECCO A/ALL 1000PKG
DETAILED DESCRIPTION	10 boxes containing 100 tubes with aluminum applicator and orange cap in bulk

PRODUCT COMPONENTS

APPLICATOR	ALUMINIUM APPLICATOR	DIMENSIONAL SPECIFICATIONS
		STICK LENGTH (*purchased product, no tolerances): 149.0 mm STICK DIAMETER HANDLE PART: 0.9 mm \pm 0.1 mm TIP DIAMETER: 3.0 mm \pm 0.5 mm FIBRE TIP LENGTH: 15.0 mm \pm 2.0 mm TOTAL LENGTH: 150.0 mm \pm 1.0 mm
CAP	ORANGE CAP	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 12.20 mm \pm 0.10 mm CAP LENGTH: 32.50 mm \pm 0.20 mm INSERTION DIAMETER: 9.30 mm \pm 0.05 mm
TRANSPORT TUBE	TUBE	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 155.00 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.5 mm EXTERNAL TUBE DIAMETER: 12.0 mm \pm 0.5 mm

This bulletin (Rev.00-05/09/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.

**TECHNICAL DATA SHEET****PACKAGING SPECIFICATION**

INNER BOX QTY	100
OUTER BOX QTY	1000
PACKAGING	BULK

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.00-05/09/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.