

Retort Stand Base

steel
thickness 8 mm
powdercoated
M 10 threaded
rubber feet

Article no. length mm width mm

2000	180	100	1100 g
2002	210	130	1700 g
2003	250	160	2500 g
2004	300	150	2800 g
2005	315	200	4000 g



© 2015

Carl Friedrich Usbeck KG
Industriestr. 12-14
42477 Radevormwald
Germany

Phone: +49 (02195) 91 18-0
Telefax: +49 (02195) 91 18-40

<https://www.usbeck.eu>