

## Base, flexible

zinc diecasting  
and brass,  
nickel plated  
vinyl coated jaws  
length 300 mm  
safety base  
Ø 82 mm

Article no.		capacity mm	
<b>2670</b>	three-prong	0-40	445 g
<b>2672</b>	round jaws	12-25	425 g



© 2015

Carl Friedrich Usbeck KG  
Industriestr. 12-14  
42477 Radevormwald  
Germany

Phone: +49 (02195) 91 18-0  
Telefax: +49 (02195) 91 18-40

<https://www.usbeck.eu>