

## Clamp, flexible with Bosshead swivel type

zinc diecasting  
and brass,  
nickel plated  
vinyl coated jaws  
length 300 mm



Article no.		capacity mm	
<b>2680</b>	three-prong	0-40	330 g
<b>2682</b>	round jaws	12-25	310 g

© 2015

Carl Friedrich Usbeck KG  
Industriestr. 12-14  
42477 Radevormwald  
Germany

Phone: +49 (02195) 91 18-0  
Telefax: +49 (02195) 91 18-40

<https://www.usbeck.eu>