

Tripod Stand

Zinc plating of tripod stands (art-no. 5131-5136) serves as temporary corrosion protection. May peel off after direct contact with open burner flame.

Choose stainless steel version (art-no. 5141-5146), alternatively.

Article no.

zinc plated steel	inside Ø mm	height mm	
5131	100	180	290 g
5132	120	210	340 g
5133	140	220	385 g
5134	160	230	400 g
5135	120	250	380 g
5136	120	260	390 g

stainless steel	inside Ø mm	height mm	
5141	100	180	290 g
5142	120	210	340 g
5143	140	220	385 g
5144	160	230	400 g
5145	120	250	380 g
5146	120	260	390 g



© 2015

Carl Friedrich Usbeck KG
Industriestr. 12-14
42477 Radevormwald
Germany

Phone: +49 (02195) 91 18-0
Telefax: +49 (02195) 91 18-40

<https://www.usbeck.eu>