




Description	10ml Tips A Type	10ml Tips B Type
Maximum Volume	10000μl	10000μl
Graduation	NA	1000μl, 2000μl, 3000μl, 4000μl, 5000μl 6000μl, 7000μl, 8000μl, 9000μl, 10000μl
Material	Virgin Polypropylene	Virgin Polypropylene
Autoclavable	Yes	Yes
Certified RNase, DNase, DNA and Pyrogen free	Yes	Yes

Hazards Data	<p>Health 0 Non-hazardous Flammability 1 Flammable (Class A fire) Reactivity 0 Non-reactive *Consult the polymer MSDS for further information on potential hazards associated with this material.</p>	
--------------	--	--

Production Standards & Procedures	<p>Product is Produced under controlled environment. All personnel are capped, masked, gloved and gowned in the assembly area. Lot size is variable due to continuous molding process. This Product is manufactured to the following specifications: Magnified visual examination for flash (+0.02in).</p>	
-----------------------------------	--	--

Packaging	<p>AT-10000A</p>  <p>Bulk Bulk-Sterile</p>	<p>AT-10000B</p>  <p>Bulk Bulk-Sterile</p>
-----------	---	---

Technical Drawing	<p>Tips 10000 A</p> 	<p>Tips 10000 B</p> 
-------------------	--	--

Order in:- AT-10000A, AT-10000B				
Volume	Bulk	Sterile	Individually Wrapper	
10000μl	B	S	W	

B=Bulk, S=Sterile, W=Wrapped