





Description	10XN Tips	10XN Filter Tips
Maximum Volume	10 $\mu$ l	10 $\mu$ l
Graduation	2 $\mu$ l, 10 $\mu$ l	2 $\mu$ l, 10 $\mu$ l
Material	Virgin Polypropylene	Virgin Polypropylene
Autoclavable	Yes	No
Certified RNase, DNase, DNA and Pyrogen free	Yes	Yes
Hazards Data	Health 0 Non-hazardous Flammability 1 Flammable (Class A fire) Reactivity 0 Non-reactive *Consult the polymer MSDS for further information on potential hazards associated with this material.	

Production Standards & Procedures	<p>Product is Produced under controlled environment.</p> <p>All personnel are capped, masked, gloved and gowned in the assembly area.</p> <p>Lot size is variable due to continuous molding process.</p> <p>Tip Filters are specifically designed to help prevent liquid bypass ,cross contamination and airborne contamination from being drawn in to the Tips, thus protecting the pipette body from contamination.</p> <p>This Product is manufactured to the following specifications:</p> <p>Magnified visual examination for flash (+0.02in).</p>
-----------------------------------	---

Low Retention Tips & Regular Tips	<p><b>Low Retention Tips</b></p> <p>Minimum retention = Maximum accuracy</p> <p>No Silicon additives = No leaching</p> <p>Manufactured from specially developed resins combined with high quality diamond finish mould creates super hydrophobic internal surface and ensures maximum recovery of sample yield.</p>	
-----------------------------------	---	--

Packaging	    
	Bulk Bulk-Sterile Racked Racked-Sterile Reloads Reloads-Sterile Reloads - Filter Individually Wrapped Tip-Sterile

Technical Drawing		
-------------------	---	---

Order in:- 10XN-S-L-F-R							<div>Tips for Perfect Pipetting</div> 
Volume			Sterile	Low Retention	Filter	Packing	
10µl	10XL	20µl				B	
200µl	200XL	50µl				R	
1000µl	1000XL	100µl				K	
10XN	20µl					W	
i.e 10µl Extra Narrow Sterile Low Retention Filter Rack Tips							
S=Sterile, L=Low Retention, F=Filter, B=Bulk , R=Rack, K=Reload, W=Wrapped							
www.accumaximum.com							