

ACCUMAX



ROBUST RESISTANT & RELIABLE

Accumax Reagent Bottle

Accumax Reagent Bottle

Reagent bottles have been developed for a wide variety of applications across many different industries. They are designed to protect the contents from contamination and evaporation, as well as to provide safe storage and easy handling.

Various types of polymers are available for reagent bottles, different polymers offer a variety of properties each having a different level of importance.

Available in various types such as wide mouth, narrow mouth, wide mouth with handle, square bottles and heavy duty bottle



Product Features



Excellent chemical resistance for a variety of lab application

Amber colour for excellent light shielding for storage of light sensitive substance



Leak proof for ultimate protection and Autoclavable*

Produced in 100K Class Clean room Facility



USFDA approved raw material of USP class VI with low leachable and extractable

Translucent Visibility and Amber colour option available



Reagent bottles with Handle for ease of carrying and Transportaion

Available in PP, HDPE and LDPE



Threads on bottle & closures have continuous straight-shoulder semibutress threads

Wide range of volume available from 4ml to 4Ltr



*All caps are made up of PP material

Ordering Information

Narrow Mouth Polypropylene Reagent Bottle - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ANB-PP-4	Reagent Bottle with cap, Narrow mouth, Polypropylene, 4ml	12	72
ANB-PP-8	Reagent Bottle with cap, Narrow mouth, Polypropylene, 8ml	12	72
ANB-PP-15	Reagent Bottle with cap, Narrow mouth, Polypropylene, 15ml	12	72
ANB-PP-30	Reagent Bottle with cap, Narrow mouth, Polypropylene, 30ml	12	72
ANB-PP-60	Reagent Bottle with cap, Narrow mouth, Polypropylene, 60ml	12	72
ANB-PP-125	Reagent Bottle with cap, Narrow mouth, Polypropylene, 125ml	12	72
ANB-PP-250	Reagent Bottle with cap, Narrow mouth, Polypropylene, 250ml	12	72
ANB-PP-500	Reagent Bottle with cap, Narrow mouth, Polypropylene, 500ml	12	48
ANB-PP-1000	Reagent Bottle with cap, Narrow mouth, Polypropylene, 1000ml	12	24
ANB-PP-2000	Reagent Bottle with cap, Narrow mouth, Polypropylene, 2000ml	6	6
ANB-PP-4000	Reagent Bottle with cap, Narrow mouth, Polypropylene, 4000ml	6	6
Narrow Mouth Polypropylene Reagent Bottle - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ANB-PP-4-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 4ml, Sterile	12	72
ANB-PP-8-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 8ml, Sterile	12	72
ANB-PP-15-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 15ml, Sterile	12	72
ANB-PP-30-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 30 ml, Sterile	12	72
ANB-PP-60-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 60ml, Sterile	12	72
ANB-PP-125-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 125ml, Sterile	12	72
ANB-PP-250-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 250ml, Sterile	12	72
ANB-PP-500-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 500 ml, Sterile	12	48
ANB-PP-1000-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 1000 ml, Sterile	12	24
ANB-PP-2000-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 2000 ml, Sterile	6	6
ANB-PP-4000-S	Reagent Bottle with cap, Narrow mouth, Polypropylene, 4000 ml, Sterile	6	6
Narrow Mouth HDPE Reagent Bottle - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ANB-HD-4	Reagent Bottle with cap, Narrow mouth, HDPE, 4ml	12	72
ANB-HD-8	Reagent Bottle with cap, Narrow mouth, HDPE, 8ml	12	72
ANB-HD-15	Reagent Bottle with cap, Narrow mouth, HDPE, 15ml	12	72
ANB-HD-30	Reagent Bottle with cap, Narrow mouth, HDPE, 30ml	12	72
ANB-HD-60	Reagent Bottle with cap, Narrow mouth, HDPE, 60ml	12	72
ANB-HD-125	Reagent Bottle with cap, Narrow mouth, HDPE, 125ml	12	72
ANB-HD-250	Reagent Bottle with cap, Narrow mouth, HDPE,250ml	12	72
ANB-HD-500	Reagent Bottle with cap, Narrow mouth, HDPE, 500ml	12	48
ANB-HD-1000	Reagent Bottle with cap, Narrow mouth, HDPE, 1000ml	12	24
ANB-HD-2000	Reagent Bottle with cap, Narrow mouth, HDPE, 2000ml	6	6
ANB-HD-4000	Reagent Bottle with cap, Narrow mouth, HDPE, 4000ml	6	6

Ordering Information

Narrow Mouth HDPE Reagent Bottle - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ANB-HD-4-S	Reagent Bottle with cap, Narrow mouth, HDPE, 4ml, Sterile	12	72
ANB-HD-8-S	Reagent Bottle with cap, Narrow mouth, HDPE, 8ml, Sterile	12	72
ANB-HD-15-S	Reagent Bottle with cap, Narrow mouth, HDPE, 15ml, Sterile	12	72
ANB-HD-30-S	Reagent Bottle with cap, Narrow mouth, HDPE, 30ml, Sterile	12	72
ANB-HD-60-S	Reagent Bottle with cap, Narrow mouth, HDPE, 60ml, Sterile	12	72
ANB-HD-125-S	Reagent Bottle with cap, Narrow mouth, HDPE, 125ml, Sterile	12	72
ANB-HD-250-S	Reagent Bottle with cap, Narrow mouth, HDPE, 250ml, Sterile	12	72
ANB-HD-500-S	Reagent Bottle with cap, Narrow mouth, HDPE, 500ml, Sterile	12	48
ANB-HD-1000-S	Reagent Bottle with cap, Narrow mouth, HDPE, 1000ml, Sterile	12	24
ANB-HD-2000-S	Reagent Bottle with cap, Narrow mouth, HDPE, 2000ml, Sterile	6	6
ANB-HD-4000-S	Reagent Bottle with cap, Narrow mouth, HDPE, 4000ml, Sterile	6	6
Narrow Mouth HDPE Reagent Bottle Amber Color- Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ANB-HD-4-A	Reagent Bottle with cap, Narrow mouth, HDPE, 4ml, Amber	12	72
ANB-HD-8-A	Reagent Bottle with cap, Narrow mouth, HDPE, 8ml, Amber	12	72
ANB-HD-15-A	Reagent Bottle with cap, Narrow mouth, HDPE, 15ml, Amber	12	72
ANB-HD-30-A	Reagent Bottle with cap, Narrow mouth, HDPE, 30ml, Amber	12	72
ANB-HD-60-A	Reagent Bottle with cap, Narrow mouth, HDPE, 60ml, Amber	12	72
ANB-HD-125-A	Reagent Bottle with cap, Narrow mouth, HDPE, 125ml, Amber	12	72
ANB-HD-250-A	Reagent Bottle with cap, Narrow mouth, HDPE, 250ml, Amber	12	72
ANB-HD-500-A	Reagent Bottle with cap, Narrow mouth, HDPE, 500ml, Amber	12	48
ANB-HD-1000-A	Reagent Bottle with cap, Narrow mouth, HDPE, 1000ml, Amber	12	24
Narrow Mouth HDPE Reagent Bottle Amber Color- Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ANB-HD-4-A-S	Reagent Bottle with cap, Narrow mouth, HDPE, 4ml, Amber, Sterile	12	72
ANB-HD-8-A-S	Reagent Bottle with cap, Narrow mouth, HDPE, 8ml, Amber, Sterile	12	72
ANB-HD-15-A-S	Reagent Bottle with cap, Narrow mouth, HDPE, 15ml, Amber, Sterile	12	72
ANB-HD-30-A-S	Reagent Bottle with cap, Narrow mouth, HDPE, 30ml, Amber, Sterile	12	72
ANB-HD-60-A-S	Reagent Bottle with cap, Narrow mouth, HDPE, 60ml, Amber, Sterile	12	72
ANB-HD-125-A-S	Reagent Bottle with cap, Narrow mouth, HDPE, 125ml, Amber, Sterile	12	72
ANB-HD-250-A-S	Reagent Bottle with cap, Narrow mouth, HDPE, 250ml, Amber, Sterile	12	72
ANB-HD-500-A-S	Reagent Bottle with cap, Narrow mouth, HDPE, 500ml, Amber, Sterile	12	48
ANB-HD-1000-A-S	Reagent Bottle with cap, Narrow mouth, HDPE, 1000ml, Amber, Sterile	12	24
Narrow Mouth LDPE Reagent Bottle - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ANB-LD-4	Reagent Bottle with cap, Narrow mouth, LDPE, 4ml	12	72
ANB-LD-8	Reagent Bottle with cap, Narrow mouth, LDPE, 8ml	12	72

Ordering Information

Narrow Mouth LDPE Reagent Bottle - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ANB-LD-15	Reagent Bottle with cap, Narrow mouth, LDPE, 15ml	12	72
ANB-LD-30	Reagent Bottle with cap, Narrow mouth, LDPE, 30ml	12	72
ANB-LD-60	Reagent Bottle with cap, Narrow mouth, LDPE, 60ml	12	72
ANB-LD-125	Reagent Bottle with cap, Narrow mouth, LDPE, 125ml	12	72
ANB-LD-250	Reagent Bottle with cap, Narrow mouth, LDPE, 250ml	12	72
ANB-LD-500	Reagent Bottle with cap, Narrow mouth, LDPE, 500ml	12	48
ANB-LD-1000	Reagent Bottle with cap, Narrow mouth, LDPE, 1000ml	12	24
ANB-LD-2000	Reagent Bottle with cap, Narrow mouth, LDPE, 2000ml	6	6
ANB-LD-4000	Reagent Bottle with cap, Narrow mouth, LDPE, 4000ml	6	6
Narrow Mouth LDPE Reagent Bottle - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ANB-LD-4-S	Reagent Bottle with cap, Narrow mouth, LDPE, 4ml, Sterile	12	72
ANB-LD-8-S	Reagent Bottle with cap, Narrow mouth, LDPE, 8ml, Sterile	12	72
ANB-LD-15-S	Reagent Bottle with cap, Narrow mouth, LDPE, 15ml, Sterile	12	72
ANB-LD-30-S	Reagent Bottle with cap, Narrow mouth, LDPE, 30ml, Sterile	12	72
ANB-LD-60-S	Reagent Bottle with cap, Narrow mouth, LDPE, 60ml, Sterile	12	72
ANB-LD-125-S	Reagent Bottle with cap, Narrow mouth, LDPE, 125ml, Sterile	12	72
ANB-LD-250-S	Reagent Bottle with cap, Narrow mouth, LDPE, 250ml, Sterile	12	72
ANB-LD-500-S	Reagent Bottle with cap, Narrow mouth, LDPE, 500ml, Sterile	12	48
ANB-LD-1000-S	Reagent Bottle with cap, Narrow mouth, LDPE, 1000ml, Sterile	12	24
ANB-LD-2000-S	Reagent Bottle with cap, Narrow mouth, LDPE, 2000ml, Sterile	6	6
ANB-LD-4000-S	Reagent Bottle with cap, Narrow mouth, LDPE, 4000ml, Sterile	6	6
Wide Mouth Polypropylene Reagent Bottle - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-PP-30	Reagent Bottle with cap, Wide mouth, Polypropylene, 30ml	12	72
AWB-PP-60	Reagent Bottle with cap, Wide mouth, Polypropylene, 60ml	12	72
AWB-PP-125	Reagent Bottle with cap, Wide mouth, Polypropylene, 125ml	12	72
AWB-PP-250	Reagent Bottle with cap, Wide mouth, Polypropylene, 250ml	12	72
AWB-PP-500	Reagent Bottle with cap, Wide mouth, Polypropylene, 500ml	12	48
AWB-PP-1000	Reagent Bottle with cap, Wide mouth, Polypropylene, 1000ml	12	24
AWB-PP-2000	Reagent Bottle with cap, Wide mouth, Polypropylene, 2000ml	6	6
AWB-PP-4000	Reagent Bottle with cap, Wide mouth, Polypropylene, 4000ml	6	6
Wide Mouth Polypropylene Reagent Bottle - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-PP-30-S	Reagent Bottle with cap, Wide mouth, Polypropylene, 30ml, Sterile	12	72
AWB-PP-60-S	Reagent Bottle with cap, Wide mouth, Polypropylene, 60ml, Sterile	12	72
AWB-PP-125-S	Reagent Bottle with cap, Wide mouth, Polypropylene, 125ml, Sterile	12	72

Ordering Information

Wide Mouth Polypropylene Reagent Bottle - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-PP-250-S	Reagent Bottle with cap, Wide mouth, Polypropylene, 250ml, Sterile	12	72
AWB-PP-500-S	Reagent Bottle with cap, Wide mouth, Polypropylene, 500ml, Sterile	12	48
AWB-PP-1000-S	Reagent Bottle with cap, Wide mouth, Polypropylene, 1000ml, Sterile	12	24
AWB-PP-2000-S	Reagent Bottle with cap, Wide mouth, Polypropylene, 2000ml, Sterile	6	6
AWB-PP-4000-S	Reagent Bottle with cap, Wide mouth, Polypropylene, 4000ml, Sterile	6	6
Wide Mouth Polypropylene Reagent Bottle - Sterile Individual Wrapped			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-PP-30-S-IW	Reagent Bottle with cap, Wide mouth, Polypropylene, 30ml, Sterile, I/W	12	72
AWB-PP-60-S-IW	Reagent Bottle with cap, Wide mouth, Polypropylene, 60ml, Sterile, I/W	12	72
AWB-PP-125-S-IW	Reagent Bottle with cap, Wide mouth, Polypropylene, 125ml, Sterile, I/W	12	72
AWB-PP-250-S-IW	Reagent Bottle with cap, Wide mouth, Polypropylene, 250ml, Sterile, I/W	12	72
AWB-PP-500-S-IW	Reagent Bottle with cap, Wide mouth, Polypropylene, 500ml, Sterile, I/W	12	48
AWB-PP-1000-S-IW	Reagent Bottle with cap, Wide mouth, Polypropylene, 1000ml, Sterile, I/W	12	24
AWB-PP-2000-S-IW	Reagent Bottle with cap, Wide mouth, Polypropylene, 2000ml, Sterile, I/W	6	6
AWB-PP-4000-S-IW	Reagent Bottle with cap, Wide mouth, Polypropylene, 4000ml, Sterile, I/W	6	6
Wide Mouth HDPE Reagent Bottle - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-HD-30	Reagent Bottle with cap, Wide mouth, HDPE, 30ml	12	72
AWB-HD-60	Reagent Bottle with cap, Wide mouth, HDPE, 60ml	12	72
AWB-HD-125	Reagent Bottle with cap, Wide mouth, HDPE, 125ml	12	72
AWB-HD-250	Reagent Bottle with cap, Wide mouth, HDPE, 250ml	12	72
AWB-HD-500	Reagent Bottle with cap, Wide mouth, HDPE, 500ml	12	48
AWB-HD-1000	Reagent Bottle with cap, Wide mouth, HDPE, 1000ml	12	24
AWB-HD-2000	Reagent Bottle with cap, Wide mouth, HDPE, 2000ml	6	6
AWB-HD-4000	Reagent Bottle with cap, Wide mouth, HDPE, 4000ml	6	6
Wide Mouth HDPE Reagent Bottle - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-HD-30-S	Reagent Bottle with cap, Wide mouth, HDPE, 30ml, Sterile	12	72
AWB-HD-60-S	Reagent Bottle with cap, Wide mouth, HDPE, 60ml, Sterile	12	72
AWB-HD-125-S	Reagent Bottle with cap, Wide mouth, HDPE, 125ml, Sterile	12	72
AWB-HD-250-S	Reagent Bottle with cap, Wide mouth, HDPE, 250ml, Sterile	12	72
AWB-HD-500-S	Reagent Bottle with cap, Wide mouth, HDPE, 500ml, Sterile	12	48
AWB-HD-1000-S	Reagent Bottle with cap, Wide mouth, HDPE, 1000ml, Sterile	12	24
AWB-HD-2000-S	Reagent Bottle with cap, Wide mouth, HDPE, 2000ml, Sterile	6	6
AWB-HD-4000-S	Reagent Bottle with cap, Wide mouth, HDPE, 4000ml, Sterile	6	6

Ordering Information

Wide Mouth HDPE Reagent Bottle - Sterile Individual Wrapped			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-HD-S-W-30	Reagent Bottle with cap, Wide mouth, HDPE, 30ml, Sterile, I/W	12	72
AWB-HD-S-W-60	Reagent Bottle with cap, Wide mouth, HDPE, 60ml, Sterile, I/W	12	72
AWB-HD-S-W-125	Reagent Bottle with cap, Wide mouth, HDPE, 125ml, Sterile, I/W	12	72
AWB-HD-S-W-250	Reagent Bottle with cap, Wide mouth, HDPE, 250ml, Sterile, I/W	12	72
AWB-HD-S-W-500	Reagent Bottle with cap, Wide mouth, HDPE, 500ml, Sterile, I/W	12	48
AWB-HD-S-W-1000	Reagent Bottle with cap, Wide mouth, HDPE, 1000ml, Sterile, I/W	12	24
Wide Mouth HDPE Reagent Bottle Amber Color - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-HD-30-A	Reagent Bottle with cap, Wide mouth, HDPE, 30ml, Amber	12	72
AWB-HD-60-A	Reagent Bottle with cap, Wide mouth, HDPE, 60ml, Amber	12	72
AWB-HD-125-A	Reagent Bottle with cap, Wide mouth, HDPE, 125ml, Amber	12	72
AWB-HD-250-A	Reagent Bottle with cap, Wide mouth, HDPE, 250ml, Amber	12	72
AWB-HD-500-A	Reagent Bottle with cap, Wide mouth, HDPE, 500ml, Amber	12	48
AWB-HD-1000-A	Reagent Bottle with cap, Wide mouth, HDPE, 1000ml, Amber	12	24
AWB-HD-2000-A	Reagent Bottle with cap, Wide mouth, HDPE, 2000ml, Amber	6	6
AWB-HD-4000-A	Reagent Bottle with cap, Wide mouth, HDPE, 4000ml, Amber	6	6
Wide Mouth HDPE Reagent Bottle Amber Color - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-HD-30-A-S	Reagent Bottle with cap, Wide mouth, HDPE, 30ml, Amber, Sterile	12	72
AWB-HD-60-A-S	Reagent Bottle with cap, Wide mouth, HDPE, 60ml, Amber, Sterile	12	72
AWB-HD-125-A-S	Reagent Bottle with cap, Wide mouth, HDPE, 125ml, Amber, Sterile	12	72
AWB-HD-250-A-S	Reagent Bottle with cap, Wide mouth, HDPE, 250ml, Amber, Sterile	12	72
AWB-HD-500-A-S	Reagent Bottle with cap, Wide mouth, HDPE, 500ml, Amber, Sterile	12	48
AWB-HD-1000-A-S	Reagent Bottle with cap, Wide mouth, HDPE, 1000ml, Amber, Sterile	12	24
AWB-HD-2000-A-S	Reagent Bottle with cap, Wide mouth, HDPE, 2000ml, Amber, Sterile	6	6
AWB-HD-4000-A-S	Reagent Bottle with cap, Wide mouth, HDPE, 4000ml, Amber, Sterile	6	6
Wide Mouth LDPE Reagent Bottle - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-LD-30	Reagent Bottle with cap, Wide mouth, LDPE, 30ml	12	72
AWB-LD-60	Reagent Bottle with cap, Wide mouth, LDPE, 60ml	12	72
AWB-LD-125	Reagent Bottle with cap, Wide mouth, LDPE, 125ml	12	72
AWB-LD-250	Reagent Bottle with cap, Wide mouth, LDPE, 250ml	12	72
AWB-LD-500	Reagent Bottle with cap, Wide mouth, LDPE, 500ml	12	48
AWB-LD-1000	Reagent Bottle with cap, Wide mouth, LDPE, 1000ml	12	24
AWB-LD-2000	Reagent Bottle with cap, Wide mouth, LDPE, 2000ml	6	6
AWB-LD-4000	Reagent Bottle with cap, Wide mouth, LDPE, 4000ml	6	6

Ordering Information

Wide Mouth LDPE Reagent Bottle - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AWB-LD-30-S	Reagent Bottle with cap, Wide mouth, LDPE, 30ml, Sterile	12	72
AWB-LD-60-S	Reagent Bottle with cap, Wide mouth, LDPE, 60ml, Sterile	12	72
AWB-LD-125-S	Reagent Bottle with cap, Wide mouth, LDPE, 125ml, Sterile	12	72
AWB-LD-250-S	Reagent Bottle with cap, Wide mouth, LDPE, 250ml, Sterile	12	72
AWB-LD-500-S	Reagent Bottle with cap, Wide mouth, LDPE, 500ml, Sterile	12	48
AWB-LD-1000-S	Reagent Bottle with cap, Wide mouth, LDPE, 1000ml, Sterile	12	24
AWB-LD-2000-S	Reagent Bottle with cap, Wide mouth, LDPE, 2000ml, Sterile	6	6
AWB-LD-4000-S	Reagent Bottle with cap, Wide mouth, LDPE, 4000ml, Sterile	6	6
Heavy Duty Bottles - PP - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AHD-PP-1000	Reagent Bottle 1000ml Heavy Duty, PP bottle with Cap	6	6
AHD-PP-2000	Reagent Bottle 2000ml Heavy Duty, PP bottle with Cap	6	6
AHD-PP-4000	Reagent Bottle 4000ml Heavy Duty, PP bottle with Cap	6	6
AHD-PP-5000	Reagent Bottle 5000ml Heavy Duty, PP bottle with Cap	4	4
Heavy Duty Bottles - PP - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AHD-PP-1000-S	Reagent Bottle 1000ml Heavy Duty, PP bottle with Cap, Sterile	6	6
AHD-PP-2000-S	Reagent Bottle 2000ml Heavy Duty, PP bottle with Cap, Sterile	6	6
AHD-PP-4000-S	Reagent Bottle 4000ml Heavy Duty, PP bottle with Cap, Sterile	6	6
AHD-PP-5000-S	Reagent Bottle 5000ml Heavy Duty, PP bottle with Cap, Sterile	4	4
Heavy Duty Bottles - HDPE - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AHD-HD-1000	Reagent Bottle 1000ml Heavy Duty, HDPE bottle with Cap	6	6
AHD-HD-2000	Reagent Bottle 2000ml Heavy Duty, HDPE bottle with Cap	6	6
AHD-HD-4000	Reagent Bottle 4000ml Heavy Duty, HDPE bottle with Cap	6	6
AHD-HD-5000	Reagent Bottle 5000ml Heavy Duty, HDPE bottle with Cap	4	4
Heavy Duty Bottles - HDPE - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
AHD-HD-1000-S	Reagent Bottle 1000ml Heavy Duty, HDPE bottle with Cap, Sterile	6	6
AHD-HD-2000-S	Reagent Bottle 2000ml Heavy Duty, HDPE bottle with Cap, Sterile	6	6
AHD-HD-4000-S	Reagent Bottle 4000ml Heavy Duty, HDPE bottle with Cap, Sterile	6	6
AHD-HD-5000-S	Reagent Bottle 5000ml Heavy Duty, HDPE bottle with Cap, Sterile	4	4
Wide Mouth PP Bottles with Handle - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ARH-PP-4000	Reagent Bottle Wide mouth Round 4000ml with cap, PP Bottle with handle	6	6
ARH-HD-4000	Reagent Bottle Wide mouth Round 4000ml with cap, HDPE Bottle with handle	6	6
ASH-PP-2000	Reagent Bottle Wide mouth Square 2000ml with cap, PP Bottle with handle	6	6

Ordering Information

Wide Mouth PP Bottles with Handle - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ASH-HD-2000	Reagent BottleWide mouth Square 2000ml with cap, HDPE Bottle with handle	6	6
ASH-PP-4000	Reagent BottleWide mouth Square 4000ml with cap, PP Bottle with handle	6	6
ASH-HD-4000	Reagent BottleWide mouth Square 4000ml with cap, HDPE Bottle with handle	6	6
Square Wide Mouth PP Bottles - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ASB-PP-60	Reagent Bottle Wide mouth Square 60ml with cap, PP Bottle	12	72
ASB-PP-175	Reagent Bottle Wide mouth Square 175ml with cap, PP Bottle	12	72
ASB-PP-250	Reagent Bottle Wide mouth Square 250ml with cap, PP Bottle	12	72
ASB-PP-500	Reagent Bottle Wide mouth Square 500ml with cap, PP Bottle	12	48
ASB-PP-1000	Reagent Bottle Wide mouth Square 1000ml with cap, PP Bottle	12	24
Square Wide Mouth HDPE Bottles - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ASB-HD-60-S	Reagent Bottle Wide mouth Square 60ml with cap, HDPE Bottle, Sterile	12	72
ASB-HD-175-S	Reagent Bottle Wide mouth Square 175ml with cap, HDPE Bottle, Sterile	12	72
ASB-HD-250-S	Reagent Bottle Wide mouth Square 250ml with cap, HDPE Bottle, Sterile	12	72
ASB-HD-500-S	Reagent Bottle Wide mouth Square 500ml with cap, HDPE Bottle, Sterile	12	48
ASB-HD-1000-S	Reagent Bottle Wide mouth Square 1000ml with cap, HDPE Bottle, Sterile	12	24
Square Wide Mouth HDPE Bottles Amber Color - Non Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ASB-HD-60-A	Wide mouth Square 60ml with cap, HDPE Bottle, Amber	12	72
ASB-HD-175-A	Wide mouth Square 175ml with cap, HDPE Bottle, Amber	12	72
ASB-HD-250-A	Wide mouth Square 250ml with cap, HDPE Bottle, Amber	12	72
ASB-HD-500-A	Wide mouth Square 500ml with cap, HDPE Bottle, Amber	12	48
ASB-HD-1000-A	Wide mouth Square 1000ml with cap, HDPE Bottle, Amber	12	24
Square Wide Mouth HDPE Bottles Amber Color - Sterile			
Cat. No.	Description	Qty./Pack	Qty./Case
ASB-HD-60-A-S	Wide mouth Square 60ml with cap, HDPE Bottle, Amber, Sterile	12	72
ASB-HD-175-A-S	Wide mouth Square 175ml with cap, HDPE Bottle, Amber, Sterile	12	72
ASB-HD-250-A-S	Wide mouth Square 250ml with cap, HDPE Bottle, Amber, Sterile	12	72
ASB-HD-500-A-S	Wide mouth Square 500ml with cap, HDPE Bottle, Amber, Sterile	12	48
ASB-HD-1000-A-S	Wide mouth Square 1000ml with cap, HDPE Bottle, Amber, Sterile	12	24

ΔACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaxlab.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.