

Fiche technique

pH-Electrode PINK 0330000

Gamme de pH:	pH 0...14
Gamme de temperature:	0...60°C
Prssion maxi:	6 bar
Matière du cormps:	POM/ABS
Connecteur:	S7 PINK
Electrolyte:	Polymer
Systeme de référence:	Ag/AgCl cartridge
Diaphragme:	2 pores
Resistance:	≤ 300 MΩ at 25°C
Membrane verre:	S
Immersion minimum:	15 mm
Point zéro:	0 +/- 20mV
Sensitivity:	57...59mV / pH at 25°C
Pente théorique:	50 μS/cm
Temps de réaction:	pH 4..7 < 30 s

