

Fiche technique

pH-Electrode PINK 0330001

Gamme de pH:	pH 0...14
Gamme de temperature:	0...60°C
Connecteur:	BNC
Electrolyte:	Gel
Systeme de référence:	Ag/AgCl
Diaphragme:	1 Ceramic
Resistance:	$\leq 300 \text{ M}\Omega$ at 25°C
Membrane verre	S
Immersion minimum:	15 mm
Point zéro:	0 +/- 20mV
Pente théorique:	57...59mV / pH at 25°C
Temps de réaction:	pH 4...7 < 30 s

