

Technical Data Sheet

pH-Electrode - 0330077

pH-Range:	pH 2...14
Temperature Range:	0...100°C
Max. working pressure:	6 bar
Shaft Material:	stainless steel
Connector:	S7 red
Electrolyte:	white Polymer
Reference system:	Ag/AgCl cartridge
Diaphragm:	2 pore
Resistance:	≤ 300 MΩ at 25°C
Membrane glass:	S
Minimum immersion depth:	15 mm
Zero-Point:	0 +/- 20mV
Sensitivity:	57...59mV / pH at 25°C
Minimum conductivity of sample:	50 µS/cm
Reaction time:	pH 4..7 < 30 s

