

## Technical Data Sheet

### pH-Electrode - 0330108

---

<b>pH-Range:</b>	pH 0...14
<b>Temperature Range:</b>	0...60°C
<b>Shaft Material:</b>	PC
<b>Connector:</b>	BNC
<b>Temperature:</b>	NTC30kΩ
<b>Electrolyte:</b>	Gel
<b>Reference system:</b>	Ag/AgCl cartridge
<b>Diaphragm:</b>	1 Ceramic
<b>Resistance:</b>	≤ 300 MΩ at 25°C
<b>Membrane glass:</b>	S
<b>Minimum immersion depth:</b>	15 mm
<b>Zero-Point:</b>	0 +/- 20mV
<b>Sensitivity:</b>	57...59mV / pH at 25°C
<b>Reaction time:</b>	pH 4..7 < 30 s

