

Technical Data Sheet

pH-Electrode (204G-DIN-BAN-PC25)

0330173

pH-Range:	pH 0...14
Temperature	Range: 0...60°C
Diameter:	12mm
Length:	250mm
Shaft Material:	PC
Connector:	DIN + Banana 4mm
Temperature:	NTC30k Ω
Electrolyte:	Gel
Reference system:	Ag/AgCl
Diaphragm:	1 Ceramic
Resistance:	$\leq 300 \text{ M}\Omega$ at 25°C
Membrane glass:	S
Minimum immersion depth:	15 mm
Zero-Point:	0 +/- 20mV
Sensitivity:	57...59mV / pH at 25°C
reaction time:	pH 4..7 < 30 s

