

Technical Data Sheet

pH-Electrode (204G-BNC-RCA-B)

0330174

pH-Range:	pH 0...14
Temperature Range:	0...60°C
Shaft Material:	PC
Diameter:	12mm
Length:	130mm
Connector:	BNC + Chinch (RCA)
Temperature:	NTC30k Ω
Electrolyte:	Gel
Reference system:	Ag/AgCl cartridge
Diaphragm:	1 Ceramic
Resistance:	$\leq 300 \text{ M}\Omega$ at 25°C
Membrane glass:	S
Minimum immersion	depth: 15 mm
Zero-Point:	0 +/- 20mV
Sensitivity:	57...59mV / pH at 25°C
Reaction time:	pH 4..7 < 30 s

