

	Product Specification	Specification Number XXXXXXXXXX FI € €
		Revision 05

Part No. ~~XX~~ FI € €

Description Open PharmaScoop® 500ml

Nominal Volume: 500ml

Material of Construction: 316L stainless steel



Method of Construction: Crevice free, all welds ground and polished

Surface Finish: Better than 0.5 microns Ra

Overall Length: 300mm


Body Length: 195mm

Body Diameter: 100mm

Nominal Weight: 330g

BSE/TSE statement: All polishing compounds used in the manufacture of this scoop are of vegetable origin

PharmaScoop is a registered trade mark belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	
		Sampling Systems Ltd.