

1. Part Number: 0415083

2. Photograph:



3. Product Information:

Description	Open PharmaScoop® 2000ml
Nominal Volume:	2000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	390mm
Body Length:	270mm
Body Width:	165mm
Nominal Weight:	690g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

