


 <b>SAMPLING SYSTEMS</b>	<b>Product Specification</b>	Specification Number	A316-1700
		Revision	02

Part No. A316-1700

Description PharmaScoop®

Nominal Volume: 1700ml

Material of Construction: 316L stainless steel



Method of Construction: Crevice free, all welds ground and polished

Surface Finish: Better than 0.5 microns Ra

Overall Length: 355mm


Body Length: 230mm

Body Diameter: 135mm

Nominal Weight: 700g

BSE/TSE statement: All polishing compounds used in the manufacture of this scoop are of vegetable origin

PharmaScoop is a registered trade mark belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 03-Dec-12
----------------------------	--	--