




Broad Spectrum Microtabs™ II

REF

953 Bronopol CAS #: 52-51-7 < 44%
Natamycin CAS #: 7681-93-8 < 2%

800 Tablets



9 5 3

Danger



Instructions for use:

1. Use one tablet per 40 mL of milk sample. Do not over-agitate. Tablet should be completely dissolved before testing.
2. Invert sample to mix.
3. In the event of a jam, invert the dispenser.
4. Store in a cool, dry place between 15 °C-25 °C (59 °F-77 °F).

Hazard Statements:

- May form combustible dust concentrations in air
- Harmful if swallowed
- Causes skin irritation
- Causes serious eye damage
- Very toxic to aquatic life

First Aid:

- **If swallowed:** Call a poison center or doctor if you feel unwell.
- **If on skin:** Wash with plenty of water.
- **If in eyes:**
 1. Rinse cautiously with water for several minutes.
 2. Remove contact lenses if present and easy to do.
 3. Continue rinsing. Immediately call a poison center or doctor.
 4. Rinse mouth.

Precautionary Statements:

- Wash hands thoroughly after handling.
- Do not eat, drink, or smoke when using this product.
- Avoid release to the environment.
- Wear eye protection and protective gloves.
- **If skin irritation occurs:** Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

ADVANCED
INSTRUMENTS

Advanced Instruments, LLC
Two Technology Way
Norwood, Massachusetts 02062, USA
781-320-9000 800-225-4034
aicompanies.com

LOT

Lot No.



Exp. Date

82PI13 Rev. 5

82PI13AW Rev0