

**TECHNICAL DATA SHEET****PRODUCT DESCRIPTION**

PRODUCT CODE	8175CS20H
GENERAL DESCRIPTION	LOOPS 1UL HARD 4000PKG
DETAILED DESCRIPTION	1 outer box contain 4 inner boxes each composed of 50 paper-film pouches each containing 20 1uL hard loops.

PRODUCT COMPONENTS

LOOP	ANSA 1 ul	SPECIFICATIONS
		TOTAL LENGTH: 199.00 mm $\pm$ 0.50 mm EXTERNAL DIAMETER RING: 3.00 mm $\pm$ 0.10 mm INTERNAL DIAMETER RING: 0.80 mm $\pm$ 0.10 mm THICKNESS RING: 1.10 mm $\pm$ 0.10 mm MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	1000
OUTER BOX QTY	4000
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	IONIZING IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.02-21/03/2023) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.