



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	179CS20
GENERAL DESCRIPTION	ANSE 10UL SOFT B/20PZ 4000PKG
DETAILED DESCRIPTION	1 outer box contain 4 inner boxes each composed of 50 paper-film pouches each containing 20 light blue soft 10uL loops.

PRODUCT COMPONENTS

LOOP	LOOP 10UL SOFT LIGHT BLUE	SPECIFICATIONS
		TOTAL LENGTH: 201.00 mm $\pm$ 1.00 mm EXTERNAL DIAMETER RING: 6.50 mm $\pm$ 0.20 mm INTERNAL DIAMETER RING: 3.90 mm $\pm$ 0.20 mm THICKNESS RING: 1.20 mm $\pm$ 0.20 mm MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	1000
OUTER BOX QTY	4000
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.02-19/05/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.