

Specification

Selective culture medium for the enumeration of coagulase-positive staphylococci from a variety of samples, acc. to standard ISO 6888-2

Presentation

10 Prepared bottles
Bottle 125 ml
with: 90 ± 3 ml

Packaging Details

1 box with 10 bottles 125 ml. Injectable cap: Plastic screw inner cap. The use of syringes needles with a diameter greater than 0.8 mm is not recommended.

Shelf Life

12 months

Storage

8-25 °C

Composition

Composition (g/l):

Casein Peptone.....	10.0
Sodium pyruvate.....	10.0
Glycine.....	12.0
Meat extract.....	5.00
Lithium chloride.....	5.00
Yeast extract.....	1.00
Agar	15.0

Description /Technique

Collect, dilute and prepare samples and volumes as required according to specifications, directives, official standard regulations and/or expected results.

Melt the medium contained in the bottles in a water bath (100 °C) or in a microwave oven, avoiding overheating, before pouring into Petri dishes when cooled to 40-45°C temperature after having added a sterile reconstituted RPF supplement vial (1 vial / 90 ml Medium Base). Do not use this base with Egg Yolk Tellurite supplement.

Once solidified on a flat surface, Spread the plates by streaking methodology or by spiral method.

Incubate the plates right side up aerobically at 37±1°C for 24-48h.

After incubation, enumerate all the colonies that have appeared onto the surface of the agar. *S.aureus* form brilliant black colonies (tellurite reduction) always surrounded by a halo of fibrin.

Each laboratory must evaluate the results according to their specifications. Presumptive isolation of *S.aureus* must be confirmed by further microbiological and biochemical tests.

Calculate total microbial count per ml of sample by multiplying the average number of colonies per plate by the inverse dilution factor if streaked a diluted sample. Report results as Colony Forming Unit (CFU's) per ml or g along with incubation time and temperature.

Quality control

Physical/Chemical control

Color : yellow

pH: 7.2 ± 0.2 at 25°C

Microbiological control

Before addition of 1 reconstituted vial of RPF supplement

Inoculate: Practical range 100 ± 20 CFU. min. 50 CFU (productivity)/ 10⁴-10⁶ (selectivity) and <100 UFC (specificity-PhEur) and ≥10³ UFC (specificity-ISO).

Analytical methodology according to ISO 11133:2014/A1:2018; A2:2020.

Aerobiosis. Incubation at 37 °C ± 1, reading after 24-48 ± 2h

Microbiological control according to ISO 11133:2014/A1:2018.

Microorganism

Escherichia coli ATCC® 8739, WDCM 00012

Stph. aureus ATCC® 25923, WDCM 00034

Stph. saprophyticus ATCC® 15305, WDCM 00159

Stph. epidermidis ATCC® 12228, WDCM 00036

Sterility control

Incubation 48 h at 30-35 °C and 48 h at 20-25 °C: NO GROWTH.

Check at 7 days after incubation in same conditions.

Growth

Inhibited

Good ≥ 50%. Black/grey colonies with opacity halo

Black/grey colonies without opacity halo

Black/grey colonies without opacity halo

Bibliography

- ATLAS R.M. & L.C. PARKS (1993) Handbook of Microbiological Media. CRC Press. London.
- BAIRD-PARKER, A.C. (1962) An improved diagnostic and selective medium for isolating coagulase-positive staphylococci. J. Appl. Bact. 25:12.
- COLIPA (1997) Guidelines on Microbial Quality Management (MQM). Brussels.
- DOWNES, F.P. & K. ITO (2001) Compendium of Methods for the Microbiological Examination of Foods. 4th ed. APHA. Washington. USA.
- EUROPEAN PHARMACOPOEIA (2007) 5thed. Suppl. 5.6 § 2.6.13 Microbiological examination of non-sterile products. EDQM. Council of Europe. Strasbourg.
- FIL-IDF 60:2001 Standard. Lait et produits à base de lait - Detection des staphylocoques à coagulase positive - Technique du nombre le plus probable. Brussels.
- ISO 5944:2001 Standard. Milk and Milk based products - Detection of coagulase positive staphylococci - MPN Technique. Geneva.
- ISO 6888-1:2021 / Adm 1:2023. Standard. Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species)- Part 1 Technique using Baird-Parker Agar medium. Adment 2: Inclusion of an alternative confirmation test using RPFA stab method.
- ISO 6888-2:2021 Standard. Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of coagulase-positive staphylococci - Part 1 Technique using rabbit plasma fibrinogen agar medium. Geneva.
- ISO 11133:2014/ Adm 1:2018. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- USP 31 - NF 26 (2008) <61> Microbial Limit Tests. US Pharmacopoeial Conv. Inc. Rockville. MD. USA.
- ZANGERL, P. & H. ASPERGER (2003) Media used in the detection and enumeration of *Staphylococcus aureus*. In Handbook.